



# Chemical resistance guide

<b>LEVEL 0</b>	<b>LEVEL 1</b>	<b>LEVEL 2</b>	<b>LEVEL 3</b>	<b>LEVEL 4</b>	<b>LEVEL 5</b>	<b>LEVEL 6</b>
< 10 min	10 > 29 min	30 > 59 min	60 > 119 min	120 > 239 min	240 > 479 min	> 480 min

## SMARTLine™ Sterile Light Nitrile 300 DI



- Length: 300 mm/ 11.8"
- Chemical performance: Type B
- Biological risk: AQL 1.0 / Level 2
- Particles level: < 3,000 particles/cm<sup>2</sup> > 0.5 µm / 2,300 particles
- Colour: White
- ESD - Static dissipative / Silicone-free
- Mechanical risk: N/A
- Applications: Cleanroom

7722-84-1 Hydrogen peroxide 30%	<b>LEVEL 6</b> 480 min DR 13%
1310-73-2 Sodium Hydroxide 40%	<b>LEVEL 6</b> 480 min DR -34%
50-00-0 Formaldehyde 37%	<b>LEVEL 6</b> 480 min DR -7%

**DISCLAIMER:** The data provided was based on gloves tested under laboratory conditions, in accordance with EN 16523-1:2015 (formerly EN 374-3:2003) and EN 374-4:2013. The information is for guidance only and may not reflect the user's application. A risk assessment should always be made by purchaser to assess the suitability of gloves for a specific application.