




Chemical resistance guide

LEVEL 0	LEVEL 1	LEVEL 2	LEVEL 3	LEVEL 4	LEVEL 5	LEVEL 6
< 10 min	10 > 29 min	30 > 59 min	60 > 119 min	120 > 239 min	240 > 479 min	> 480 min

SMARTLine™ Sterile Latex 600 DI



- Length: 600 mm/ 23.6"
- Chemical performance: Type B
- Biological risk: AQL 1.0 / Level 2
- Particles level: < 3,000 particles/cm² > 0.5 µm / 2,300 particles
- Colour: Natural colour
- Silicone-free
- Mechanical risk: N/A
- Applications: Cleanroom

7722-84-1 Hydrogen peroxide 30%	LEVEL 6 480 min DR -44%
1310-73-2 Sodium Hydroxide 40%	LEVEL 6 480 min DR -41%
50-00-0 Formaldehyde 37%	LEVEL 6 480 min DR -46%

DISCLAIMER: The data provided was based on gloves tested under laboratory conditions, in accordance with EN 16523-1:2015 (formerly EN 374-3:2003) and EN 374-4:2013. The information is for guidance only and may not reflect the user's application. A risk assessment should always be made by purchaser to assess the suitability of gloves for a specific application.