

pipetman®

THE PIPETTING STANDARD

**RECOGNIZED RELIABILITY AND PROVEN DURABILITY COMBINED
WITH EASY PIPETTING AND HANDLING**



Legendary Performance

Accurate and precise, PIPETMAN® air-displacement pipettes are recognized as one of the world standard. They come in a large range of volumes, including single and multichannel models, allowing you to handle various sample types with confidence.

Easy Pipetting

PIPETMAN feels like PIPETMAN® Classic but with lower forces for enhanced pipetting habits with a standard pipette.

Built-to-Last and Convenient

Made of PVDF and stainless steel, PIPETMAN offers premium durability. Built-to-last, it is the economical solution for years of pipetting.

Simple to Maintain

Routine cleaning and maintenance are all that is required to keep PIPETMAN in top condition for years. The PIPETMAN design allows for easy access to user-serviceable parts, along with a fully autoclavable tip holder and ejector.



Choose Your Pipette Stand

Different pipette stands and holders allow you to store pipettes vertically, protecting them from contamination.

CARROUSEL™



TRIO™



Universal
Multichannel
Stand



SINGLE™




Pipette Stands

Description	Part Number
CARROUSEL™ Pipette Stand (7 pipettes)	F161401
TRIO™ Pipette Stand (3 pipettes)	F161405
Universal Multichannel Stand	F161417
SINGLE™ Pipette Holder	F161406

Choose Your Pipette Kits

PIPETMAN single channel models are available in kits, offering the essential pipetting tools to start an experiment right away.



PIPETMAN® Kits

Model	Part Number	Models
PIPETMAN® 4-Pipette Kit	F167360	P2 P20 P200 P1000
PIPETMAN® Starter Kit	F167900	P20 P200 P1000
PIPETMAN® Micro-volume Kit	F167800	P2 P10 P100
PIPETMAN® Large-volume Kit	F167920	P5000 P10mL

PIPETMAN® Single Channel

Model	PIPETMAN® DIAMOND Tips	Part Number	Volume
P2	D10 DL10	DF10ST DFL10ST	F144054M 0.2-2 µL
P10	D10 DL10	DF10ST DFL10ST	F144055M 1-10 µL
P20	D200	DF30ST	F144056M 2-20 µL
P100	D200	DF100ST	F144057M 10-100 µL
P200	D200 D300	DF200ST DF300ST	F144058M 20-200 µL
P1000	D1000 D1200	DF1000ST DF1200ST	F144059M 100-1000 µL
P5000*	D5000		F144066 500-5000 µL
P10mL*	D10mL		F144067 1-10 mL


PIPETMAN® Multichannel

P8x10 P12x10	D10 DL10	DF10ST DFL10ST	F144068 F144069	1-10 µL
P8x20 P12x20	DL10 D200	DFL10ST DF30ST	F144070 F144071	2-20 µL
P8x200 P12x200	D200 D300	DF200ST DF300ST	F144072 F144073	20-200 µL
P8x300 P12x300	D200 D300	DF200ST DF300ST	F144074 F144075	30-300 µL

Models equipped with a tip ejector are only sold with a stainless steel tip ejector. Plastic tip ejector is available as an accessory. Refer to *PIPETMAN® User's Guide LT801122*.

*P5000 and P10mL do not have tip ejectors.

Gilson maximum permissible errors are guaranteed only when PIPETMAN® pipettes are used with the recommended PIPETMAN DIAMOND Tips.



PIPETMAN DIAMOND Tips are designed to fit your PIPETMAN perfectly and guarantee maximum accuracy and precision.

For ordering information, visit www.gilson.com

gilson.com/contactus

LT800591/F | ©2023 Gilson S.A.S.
All rights reserved. All trademarks
are the property of Gilson, Inc.
and its subsidiaries.

