

BRAND GMBH + CO KG | Otto-Schott-Str. 25 | 97877 Wertheim | Germany

BRAND GMBH + CO KG  
P.O. Box 11 55  
97861 Wertheim  
Germany  
T +49 9342 808 0  
F +49 9342 808 98000  
info@brand.de  
www.brand.de

Wertheim, 28/08/2023

## Product Information: Dispensette®, applicability of CE marking

Dear valued customer and business partners,

We hereby inform you that the bottle-top dispensers Dispensette® S, Dispensette® S Organic and Dispensette® S Trace Analysis are not subject to any EC directive which comprises a CE marking. It is therefore neither required nor permitted to apply a CE marking on these bottle-top dispensers.

We confirm that the bottle-top dispensers Dispensette® S, Dispensette® S Organic and Dispensette® S Trace Analysis are produced within the specifications and error limits of DIN EN ISO 8655 and are DE-M marked according to the German Measurement and Calibration Regulation. Further details are given in the current manual provided on our webpage [www.brand.de](http://www.brand.de).

With best regards,

BRAND GMBH + CO KG

i.V.



Dr. Antonio Romaguera

- Head of Product Management-

i.A.



Dr. Patrick Jost

- Product Manager Industrial Products -