


Schraubkappe GL45/ Cap GL45
 PP blue: 0668112-GZA-0GL45-..
 Ausgießring / pouring ring
 PP blue: 0672101-GZA-0GL45-..

Gewinde/Screw-thread
 DIN 168-GL 45x4; 3 turns
 0303803-GZA-01694-..

Rohkörper/blank
 D503301-GZA-04171-..

L	Add English text	10-JUL-2023	Weindorf
K	Warnhinweis hinzu	02-MAR-18	Becker
J	Logotausch DWK	24-JAN-18	Becker
Index	Aenderung / changes	Datum / date	Name



Use blank 0311009-GZA-04171-..

Dimensions not fitted with specific tolerances are to be seen as geometric dimensions only and are not subject to any tests/checks during production. These dimensions are approximate values only and could vary in the final products.

The reproduction, distribution and utilization of this document as well as the communication of its contents to others without express authorization is prohibited. Offenders will be held liable for the payment of damages. All rights reserved in the event of the grant of a patent, utility model or design.

zul. Abweichung allowed dev.		Werkstoff / material DURAN		Stück quantity	Gewicht weight	Pos. position
Oberfläche surface		Benennung / title Laboratory Bottle 20,000 mL, GL45				Maßstab scale 1 : 2
Datum / date	Name	Art.-No. see tab.				
24.01.2018	Becker					
bearb. drawn		Form-Nr. mould number		4171/20L		A3
geprüft checked		Zeichnungsnummer / drawing number		0311006-GZA-04171		Revision Blatt/sheet L 1/1
		Ers. / repl.f.		Ers. / repl.d.		
		Schutzvermerk ISO 16016 beachten / Refer to protection notice ISO 16016				INV

21 806 91 04	Flasche braun ohne Kappe und Ring / amber bottle without cap and ring
21 801 91 57	mit blaue PP Kappe und Ausgießring / blue PP cap and pouring ring
21 801 91 08	Flasche klar ohne Kappe und Ring/ clear bottle without cap and ring
Art.-No.	Benennung/title