

DWK Life Sciences GmbH

Hattenbergstraße 10
D-55122 Mainz
Germany

Phone +49 (0)6131-1445 4131
Fax +49 (0)6131-1445 4016
www.DWK-LifeSciences.com

Date: 05.11.2025

To whom it may concern,

DWK Order code	Description
218019251	DURAN® Original GL 45 Laboratory Bottle, clear, with screw cap and pouring ring, PP, blue, 25,000 mL

1. Recommended Shelf-life

We suggest a 5-year shelf life from the date of manufacture (indicated by 8-digit Retrace code (Glass Bottles) if stored under the recommended storage conditions.

2. Recommended storage conditions for DURAN Labware Products

According to USP <659> "Packaging and Storage Requirements", articles must be protected from moisture, freezing, excessive heat and from light during shipping and distribution.

Please do not hesitate to contact me if you require any further information.

Yours sincerely,



i.V Alistair Rees
Product Manager Labware

DWK Life Sciences GmbH

Hattenbergstraße 10
D-55122 Mainz
Germany

Phone +49 (0)6131-1445 4131
Fax +49 (0)6131-1445 4016
www.DWK-LifeSciences.com

Date: 05.11.2025

To whom it may concern,

Product Manufacturing Location

Order code	Description
218019251	DURAN® laboratory bottle clear with GL 45 thread complete with Polypropylene pouring ring and Polypropylene screw cap 25,000 ml

We hereby confirm that the article listed above is made from DURAN® borosilicate 3.3 glass and that:

- The supplier is DWK Life Sciences GmbH
- The country of origin is Croatia
- And the processing steps and the addresses of the manufacturing sites are:

Process step	Process step description	Manufacturing Location	Ownership
1	Glass melting tank	DURAN d.d. Mazuraniceva 3, 52100 Pula, Croatia.	DWK Life Sciences
2	Blow moulding of DURAN® glass laboratory bottle body.	DURAN d.d. Mazuraniceva 3, 52100 Pula, Croatia.	DWK Life Sciences
3	The attachment of the GL 45 glass thread to main bottle body.	DURAN d.d. Mazuraniceva 3, 52100 Pula, Croatia.	DWK Life Sciences
4	Polypropylene pouring ring and Polypropylene screw cap	Approved DWK Life Sciences GmbH supplier, Thuringia, Germany	3 rd party supplier
5	Final assembly of components and packing of product. Final quality release of product. CoC issue.	DURAN d.d. Mazuraniceva 3, 52100 Pula, Croatia.	DWK Life Sciences

Please do not hesitate to contact me if you require any further information.

Yours sincerely,



i.V. Alistair Rees Product Manager Labware

DWK Life Sciences GmbH

Managing Directors: Armin Reiche, Michael Mora, Bernhard Scherer

Headquarters: Wertheim • District Court Mannheim • HRB 570 185 • VAT No. DE811119982