



Thread
DIN 168-GL 45x4, 3 turn
dwg No. 0303803-GZA-01694

Capacity mark

225,5 ± 1,5

~2 vaulted



109 76 61	brown, light blue cap + transp. ring	1519/10D
109 76 62	brown, red cap + red ring	1519/10Red
22 806 54 07	brown, without cap + ring	1519/10
22 801 54 51	light blue cap + transp. ring	1516/10D
22 801 54 19	red cap + red ring	1516/10Red
22 801 54 02	without cap + ring	1517/10
20 801 54 42	blank bottle with print	N/A
Art.-No.	title	Code

Weitergabe sowie Vervielfältigung dieses Dokuments, Verwertung und Mitteilung seines Inhalts sind verboten, soweit nicht ausdrücklich gestattet. Zuwiderhandlungen verpflichten zu Schadenersatz. Alle Rechte für den Fall der Patent-, Gebrauchsmuster- oder Geschmacksmustereintragung vorbehalten.

The reproduction, distribution and utilization of this document as well as the communication of its contents to others without express authorization is prohibited. Offenders will be held liable for the payment of damages. All rights reserved in the event of the grant of a patent, utility model or design.

K	was www.scilabware.com	29-SEP-2020	Weindorf
I	Silkscreen changed was 0140101-SZ1-C	01-DEC-2014	Weindorf
Index	Aenderung / changes	Datum / date	Name
zul. Abweichung allowed dev.	Werkstoff / material	Stück quantity	Gewicht weight
Oberfläche surface	Benennung / title	Pos. position	
Datum / date	Name	Laborflasche mit Gewinde; 1000ml	
bearb. drawn	01.12.2014	Weindorf	DWK Life Sciences Limited
geprüft checked			Art.-No. see Tab.
Form-Nr. mould number		1694/.1L	A3
Zeichnungsnummer / drawing number		0303904-GZA-01694	Revision Blatt/sheet
Ers. / repl.f.		K	1/1
Schutzvermerk ISO 16016 beachten. / Refer to protection notice ISO 16016.		INV	

DURAN WHEATON KIMBLE
DURAN Produktions GmbH & Co. KG
www.DWK-LifeSciences.com