

# ISO17025 Pre-Calibrated Pipettors

## Corning® Lambda™ Plus and Axygen® Axypet®



### Product Information:

- ISO17025 Pre-Calibrated Corning Lambda Plus and Axygen Axypet (single-, 8- and 12-channel) – new set of SKU's delivered with ISO17025 calibration certificate recognized globally
- Choice of the testing methods (separate SKU's)
  - Standard Procedure: 10 measurements at 3 volumes (nominal volume, 50% of nominal volume, 10% of nominal volume or minimal volume)
  - Alternative Procedure: 4 measurements at 3 volumes (nominal volume, 50% of nominal volume, 10% of nominal volume or minimal volume)
- ISO17025 certificate issued by the Warsaw Plant Accredited Laboratory (PZ HTL SA). Accreditation body: PCA (Polish Centre for Accreditation) signatory to ILAC MRA which ensures recognition of calibration certificates worldwide (A2LA, IAS, NVLAP in USA)

### Applications:

- All pipetting tasks in microbiology, molecular biology, clinical assays, immunoassays, cell culture

### Prospecting:

- Is Standard Quality Control Certificate Satisfactory to you?
- Do you need ISO 17025 Accredited Calibration Certificate for your pipettors? – We can offer Pre-Calibrated pipettors (brand new pipettors supplied with ISO17025 certificate)

### Target customers:

- All laboratories that purchase pipettors: pharmaceutical, government, environmental, food and beverage
- Not all customers know about, understand or comply with ISO 17025. Look for those that adopted ISO 17025 or equivalent rigorous system

### Tools:

Sell Sheets:

- CLS – A – EQ – 133 A4 – Axygen Axypet
- CLS – EQ – 128 A4 – Corning Lambda Plus

