

Corning® Octagonal PET Bottles

For storage and sampling of solutions

CORNING

New polyethylene terephthalate (PET) bottles are ideal for storage and sampling of media, biologicals, and other aqueous solutions. The PET resin makes the bottle lightweight, break resistant and as transparent as glass.

Features

- ▶ Available in six sizes: 30 mL, 60 mL, 125 mL, 250 mL, 500 mL and 1L
- ▶ Packaged in convenient, shrink-wrapped trays within an outer bag to assure cleanliness
- ▶ Break-resistant bottles are an ergonomic and easy to handle alternative to glass
- ▶ Bottles are sterilized to 10^{-6} SAL and certified nonpyrogenic.
- ▶ Leakproof, tamper-evident HDPE screw cap
- ▶ Molded graduations for accurate measurements
- ▶ USP class VI compliant bottles are validated noncytotoxic and nonhemolytic

New!



Corning® Octagonal PET Bottles Technical Information

Visit the Technical Resources section of our website at <http://www.corning.com/lifesciences> for these and other technical resources:

- ▶ Properties of Corning Plasticware
- ▶ Chemical Compatibility of Corning Plasticware

Corning Octagonal PET Bottles Ordering Information

Cat. No.	Capacity (mL)	Bottle Material	Screw Cap Color	Neck I.D. (mm)	Screw Cap (mm)	Qty/Tray	Qty/Cs
New 431729	30	PET*	Natural	17.5	21.0	24	120
New 431730	60	PET	Natural	17.5	21.0	24	96
431731	125	PET	Natural	31.7	34.7	24	48
431732	250	PET	Natural	31.7	34.7	24	48
431733	500	PET	Natural	31.7	34.7	12	24
431734	1000	PET	Natural	31.7	34.7	12	24

*Polyethylene terephthalate

For a complete list of Corning storage bottles, search for “Storage Bottle Selection Guide” in the “Other Site Sections” of our website www.corning.com/lifesciences.



Molded graduations for accurate measurements



Tamper-evident screw cap



Bottles are packaged in convenient, shrink-wrapped trays within an outer bag to assure cleanliness.

For additional Corning product, technical, or distributor information, please e-mail us at cctech@corning.com, visit our website www.corning.com/lifesciences or call + 31 20 659 6051 or contact your local Corning sales office listed below.

Corning Incorporated Life Sciences

Corning BV
Fogostraat 12
1060 LJ Amsterdam
The Netherlands
Phone: +31 0 20 659 60 51
Fax: +31 20 659 76 73
cceurnl@corning.com
www.corning.com/lifesciences

Worldwide Support Offices

ASIA/PACIFIC

Australia/New Zealand
t 0402-794-347

China
t 86 21 2215 2888
f 86 21 6215 2988

India
t 91 124 4604000
f 91 124 4604099

Japan
t 81 3-3586 1996
f 81 3-3586 1291

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6733-6511
f 65 6861-2913

Taiwan
t 886 2-2716-0338
f 886 2-2516-7500

EUROPE

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 31 20 655 79 28
f 31 20 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

All Other European Countries

t 31 (0) 20 659 60 51
f 31 (0) 20 659 76 73

LATIN AMERICA

Brasil
t (55-11) 3089-7419
f (55-11) 3167-0700

Mexico
t (52-81) 8158-8400
f (52-81) 8313-8589

The Corning Family of Brands

CORNING **costar** PYREX® **AXYGEN**® GOSELIN™