

Corning® 1536 Well Cyclic Olefin Copolymer (COC) Microplates

CORNING


Key Features

- ▶ Cyclic Olefin Copolymer (COC) material
- ▶ Flat top for improved sealing
- ▶ 1536 Well Low Base microplates (black with clear bottom)
- ▶ Ultra clear 127 µm film thickness
- ▶ 1536 Well Hi Base solid color microplates (black or white)
- ▶ Variety of treatment options
 - Tissue culture treated
 - Corning CellBIND® Surface
 - Poly-D-Lysine
 - Nonbinding Surface (NBS™)
 - Not treated

- ▶ Bar coded
- ▶ Custom bar codes available for compatibility with the Kalypsys System

Key Benefits

- ▶ Low autofluorescence
- ▶ Broad chemical resistance including DMSO and alcohol, which allows for assay step reduction
- ▶ High mechanical stability
- ▶ Optimized for flatness and uniformity
- ▶ Low birefringence



Corning® 1536 Well Cyclic Olefin Copolymer Microplates Ordering Information

| Cat. No. | Description | Treatment | Sterile | Qty/ Pk | Qty/ Cs |
|----------|---|---------------------------|---------|------------|------------|
| 4560 | 1536 Well Microplate, black with clear bottom | Not treated | No | 20 | 100 |
| 4561 | 1536 Well Microplate, black with clear bottom | Tissue culture treated | Yes | 20 | 100 |
| 4563 | 1536 Well Microplate, black with clear bottom | Corning CellBIND® Surface | Yes | 20 | 100 |
| 4564 | 1536 Well Microplate, black with clear bottom | Poly-D-Lysine | No | 20 | 100 |
| 4565 | 1536 Well Microplate, solid black | Not treated | No | 10 | 50 |
| 4566 | 1536 Well Microplate, solid black | Tissue culture treated | Yes | 10 | 50 |
| 4567 | 1536 Well Microplate, solid black | NBS™ | No | 10 | 50 |
| 4568 | 1536 Well Microplate, solid black | Corning CellBIND Surface | Yes | 10 | 50 |
| 4570 | 1536 Well Microplate, solid white | Not treated | No | 10 | 50 |
| 4571 | 1536 Well Microplate, solid white | Tissue culture | Yes | 10 | 50 |
| 4572 | 1536 Well Microplate, solid white | NBS | No | 10 | 50 |
| 4573 | 1536 Well Microplate, solid white | Corning CellBIND Surface | Yes | 10 | 50 |



www.corning.com/lifesciences/solutions

At Corning, we continuously strive towards improving efficiencies and developing new products and technologies for life science researchers. From assay preparation to storage, our technical experts understand your challenges and your increased need for high-quality products.

It is this expertise, plus a 160-year legacy of Corning innovation and manufacturing excellence, that puts us in a unique position to be able to offer a beginning-to-end portfolio of high-quality, reliable life sciences consumables.

For additional product or technical information, and a complete listing of our International Offices and Distributors, please visit www.corning.com/lifesciences. You can contact our **Scientific Support Team** at ScientificSupportEMEA@corning.com or through your local Corning sales office - see phone numbers listed below.

Corning Incorporated Life Sciences

Corning BV
Fogostraat 12
1060 LJ Amsterdam
The Netherlands
Phone: +31 0 20 659 60 51
Fax: +31 20 659 76 73

www.corning.com/lifesciences

European Support Offices

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 31 20 655 79 28
f 31 20 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

**All Other European
Countries**
t 31 (0) 20 659 60 51
f 31 (0) 20 659 76 73

CORNING | **FALCON**® | **AXYGEN**® | **PYREX**® | **GOSSSELIN**™

For a listing of trademarks, visit us at www.corning.com/lifesciences/trademarks.
All other trademarks in this document are the property of their respective owners.