

NOT TO SCALE

Pub No. MAN0014410

127.00 ±0.25

90.00 ±0.25

18.00

17.00 ±0.10
(x48)

85.00 ±0.25

11.00 ±0.25

9.00

8.00 ±0.10
(x48)

63.00 ±0.25

56.70 ±0.25

3.50 ±0.25
39.00 ±0.25

SECTION X-X

18.50 ±0.25

52.95 ±0.10

SECTION Y-Y

TITLE

6.0ml Storage Plate

DRG. NO.

AB-0988

ISSUE

B

DATE

24/5/02

Copyright Abgene®
This drawing and/or any technical information relating to it is issued in strict confidence and may not be reproduced or used for manufacture without written permission. Upon receipt of which it is the users responsibility to seek confirmation of any dimensions and production tolerances.

All dimensions in mm