


NOT TO SCALE


GENERAL TOLERANCE $\pm 1.0\text{mm}$

Copyright Thermo Fisher Scientific
This drawing and/or any technical information relating to it is issued in strict confidence and may not be reproduced or used for manufacture without written permission. Upon receipt of which it is the users responsibility to seek confirmation of any dimensions and production tolerances.



FOIL

FOIL 0.042 THICK, ACRYLIC ADHESIVE 0.016 THICK



RELEASE LINER

RELEASE LINER 0.052 THICK

PERFORATION DETAILS

1.0mm Cut/0.50mm Tie

TITLE

Self Adhesive PCR Foil

DRG. NO.

AB-0626

ISSUE

A

DATE

1/5/15

All dimensions in mm