

This drawing and/or any technical information relating to it is issued in strict confidence and may not be reproduced or used for manufacture without written permission. Upon receipt of which it is the users responsibility to seek confirmation of any dimensions and production tolerances.

127.76 ±0.25 NOT TO SCALE

99.00 ±0.25

14.38 ±0.25

9.00

X

63.00 ±0.25


11.24 ±0.25

9.00

2 LOCATION RIBS typ

X

85.48 ±0.25



SECTION X-X

SOLID END TUBE

28.14

31.50

38.50

16.50

4 HOLES Ø6.5 thro'

All dimensions in mm

15.00

15.00

6.75

RACK

LID

TITLE

0.65ml Cluster Tube Rack

DRG. NO.

AB-1256

ISSUE

B

DATE

18/2/08