

NOT TO SCALE
127.76 ±0.25

99.00 ±0.25

14.38 ±0.25

9.00

15.50
±0.25

85.48 ±0.25

63.00 ±0.25

11.24 ±0.25

9.00

4.50 ±0.25

25.00

25.00

27.00

4.00 (x8)

25.00

25.00


ø5.50
±0.10

0.50

15.00 ±0.10

VIAL FORM 96pos

PLATES STACKED



SECTION X-X

Copyright Abgene®
This drawing and/or any technical information relating to it is issued in strict confidence and may not be reproduced or used for manufacture without written permission. Upon receipt of which it is the users responsibility to seek confirmation of any dimensions and production tolerances.

TITLE

Thermo-Fast® 96,
Skirted PCR Plate

DRG. NO.

AB-0800

ISSUE

D

DATE

7/5/02

All dimensions in mm