

Wonder Spin



One Spin

Single Tube Spin-able Centrifuge



Multi Spin

Cordless Multi-purpose Centrifuge

■ SPECIFICATIONS

Single Tube Spin-able Centrifuge

Model	One Spin
Max. Capacity	2.0 ml × 4 microtubes
Tubes	2.0 ml × 4 microtubes 1.5 ml × 4 microtubes
Speed	4,000rpm (±20%)
RCF	450 ~ 1,000G (2.0 ml) 420 ~ 950G (1.5 ml)
Safety Devices	Lid opening interlocking brake, Overcurrent fuse, Motor overcurrent detector (Motor lock, Unbalance)
Electric Isolation	ClassII
Environmental Requirements	Ambient temperature range : from 5 to 35°C, Relative humidity : 75RH or less, Altitude : not exceeding 2000 m, For indoor use only, overvoltage category : II, pollution degree : 2
Dimensions	105W×130D×96H mm
Net Weight	380 g
Power Requirements	Dedicated AC adapter Input : 1-Phase AC100-240V 50/60Hz 0.6A
Accessories Included	•AC adapter (6A-181WP05, 5 types of plugs) × 1set •Operator's manual × 1 copy

Cordless Multi-purpose Centrifuge

Model	Multi Spin	
Rotor	MA-1	PA-1
Max. Capacity	2.0ml × 8 microtubes	0.2 ml × 8 PCR tubes (2 strips)
Tubes	2.0 ml × 8 microtubes 1.5 ml × 8 microtubes	0.2 ml × 8 PCR tubes (2 strips)
Speed	6,000rpm (±20%)	6,200rpm (±20%)
RCF	1,310 ~ 2,960G (2.0 ml) 1,260 ~ 2,840G (1.5 ml)	Max. 1,200 ~ 2,690G Min. 810 ~ 1,820G
Timer	6 minutes (Fixed)	
Safety Devices	Lid opening interlocking brake, Overcurrent fuse, Motor overcurrent detector (Motor lock, Unbalance)	
Electric Isolation	ClassII (with AC adapter)	
Environmental Requirements	Ambient temperature : from 5 to 35 °C, Relative humidity : 75%RH or less, Altitude : not exceeding 2,000 m, For indoor use only, overvoltage category : II, pollution degree : 2	
Dimensions	150W×163D×113H mm	
Net Weight	600 g (excluding rotors and batteries)	
Power Requirements	Dedicated AC adapter Input : 1-Phase AC100-240V 50/60Hz 0.6A or 4×AA size alkaline battery	
Accessories Included	• AC adapter (6A-181WP05, 5 types of plugs) × 1 set • Microtubes rotor MA-1 × 1pc. • Tube-strips rotor PA-1 × 1pc. • Operator's manual × 1copy	

- STANDARD ACCESSORY
- AC adapter
 - 5 type power plug adapters



Sales Office:
TOMY DIGITAL BIOLOGY CO., LTD.
 3-14-17 Tagara, Nerima-ku, Tokyo 179-0073, Japan
 e-mail: info@digital-biology.co.jp
 URL: http://www.digital-biology.co.jp
 phone: +81-3-5971-8160 fax: +81-3-3970-6036

TOMY SEIKO CO., LTD.

Manufacturer:
TOMY KOGYO CO., LTD.
 3-14-17 Tagara, Nerima-ku, Tokyo 179-0073, Japan

All TOMY products have a limited one-year warranty.
 Specifications are subject to change according to product advancement.
 Tomy and Digital Biology is registered trademark of Tomy Seiko Co.,
 Ltd. and Tomy Digital Biology Co., Ltd. Copyright 2016,
 Tomy Seiko and its subsidiaries. Printed in Japan.

3K0528②NC

Wonder Spin

Mini Personal Centrifuges

Save the minutes from daily routine work
 at small sample centrifugation.



Spin-able with
 a single tube!
One Spin
 Single Tube Spin-able
 Centrifuge



Chargeable
 battery drive!
Multi Spin
 Cordless Multi-purpose
 Centrifuge

TOMY

Wonder Spin


One Spin

Single Tube Spin-able Centrifuge




Active in anywhere at Bio-laboratory World-wide


Two way drive, AC power drive and cordless battery drive, expand usability. Flexible mobility at biosafety cabinet & clean bench use. 5 power plug adapters included to fit to world-wide universal power line(AC100-240V). Comply with CE standard.




- STANDARD ACCESSORY
- AC adapter
- 5 type power plug adapters



Two rotors included to spin 2.0ml tube and PCR strip



Just one push action to Close/Open




Motor over current & Imbalance detect system


6 minutes Timer equipped

Quick spin-down


A single microtube is able to spin without balance load. (Patent pending)



Lid close to spin, Release hand to stop



Easy tube unloading



Fits to universal power line

5 power plug adapters included to fit to universal power line (AC100-240V). Comply with CE standard.



Spin-column or Screw capped tube available

Compact Design

Small foot print to save your bench space.

Multi Spin

Cordless Multi-purpose Centrifuge



Wonder Spin