



**Export and Production
- since 1949 -**

Heinz Herenz Medizinalbedarf GmbH • Rudorffweg 10 • 21031 Hamburg, Germany

CERTIFICATE OF CONFORMANCE

Item: Steristoppers® No. 45 P
Item No.: 1014501

Material: 2-ply cellulose made of 100% chlorine-free cellulose
Sinking time: Max. 5 sec.
Ability of water holding: Min. 13 g/g

Glue: Glue on basis of methyl-cellulose sprayed on the cellulose while twisting it

Shape: Mould twisting process

Frequency of use: For single use

Measuring tolerances: Due to softness of the material

Sterilization: Autoclavable or dry heat up to 180-200°C

Aspiration: Steristoppers® are air permeable

Experience: Steristoppers® have been used to seal tissue cultures for many years in bacteriological institutes all over the world to entire satisfaction

Patents: German Patent No. 1467937 and foreign patents registered.

HS-code: 4823 7090

Manufacturing Date: week 32 in year 2024

Manufacturing Lot No.: 32/24
printed on each carton: week/year

Shelf Life: N/A

Inga Kieper (Quality Manager)

Date: 09.08.2024

**** This is a computer generated copy of the original document. Proof on file. No signature required. ****

www.herenz.de
Phone +49-040-739-20-40
Fax +49-040-739 20-439
Email hhh@herenz.de
ILN 40 49434 00000 8

Managing Directors
Dieter Cyll, Antje Herenz-Cyll
Comm. Rg. Hamburg HRB 48884
Tax No. 44/729/00248
V.A.T. No. DE 118 665 646

Postbank Hamburg
Sort Code 200 100 20
Account No. 6846 46 200
IBAN DE97200100200684646200
Swift-Code PBNKDEFF

Deutsche Bank AG
Sort Code 200 700 00
Account No. 55 00 251
IBAN DE82200700000550025100
Swift-Code DEUTDEHH

Hamburger Sparkasse
Sort Code 200 505 50
Account No. 1085 243 002
IBAN DE87200505501085243002
Swift-Code HASPDEHH