


CATALOG NO.....:	2226			
SIZE CODE.....:	-0020		-0050	
CARBOY MAT'L.....:	PP		PP	
CLOSURE MAT'L.....:	WHITE PP		WHITE PP	
SIZE (NOM VOL L/GAL):	10/2.5		20/5	
NOM. BRIM CAP. (L)...:	12		22.5	
NECK FINISH (mm).....:	83B		83B	
THREAD (T.P.I.).....:	2		2	
NOM. WEIGHT W/ CLOSURE(gms).....:	1560		2800	
	INCHES	MM	INCHES	MM
A (NECK I.D.).....:	2.52 ±.10	64.0 ±2.5	2.52 ±.10	64.0 ±2.5
B (HT. W/ CLOSURE)...:	15.30 ±.18	388.6 ±4.6	20.7 ±.30	526 ±8
C (HT. W/O CLOSURE)..:	14.82 ±.18	376.4 ±4.6	20.2 ±.30	513 ±8
D (DIA).....:	9.85 ±.12	250.2 ±3.0	11.10 ±.20	281.9 ±5.1
W.T.(NOM WALL THK)...:	0.14	3.6	0.15	3.8



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE PB	APP'D GRP	PSS 00-2226-XXXX-1	REV. 0606
CK'D KR			

DESCRIPTION:
VACUUM ROUND CARBOY WITH HANDLES