


| CATALOG NUMBER | DESCRIPTION                        |
|----------------|------------------------------------|
| 1200-2300      | CONVOLUTED BOTTOM, PP, NON-STERILE |
| 1200-2301      | CONVOLUTED BOTTOM, PP, STERILE     |

|                                   |           |       |
|-----------------------------------|-----------|-------|
| CATALOG NUMBER.....:              | SEE ABOVE |       |
| NOM. CAP (mL).....:               | 300       |       |
| MATERIAL.....:                    | PP        |       |
|                                   | INCH      | MM    |
| A (BASE LENGTH).....:             | 5.03      | 127.8 |
| B (BASE WIDTH).....:              | 3.365     | 85.5  |
| C (TOP LENGTH).....:              | 4.88      | 123.8 |
| D (TOP WIDTH).....:               | 3.22      | 81.7  |
| E (STEP HEIGHT).....:             | .10       | 2.5   |
| F (OVERALL HEIGHT).....:          | 1.57      | 39.9  |
| G (BOTTOM OF WELL DISTANCE).....: | .10       | 2.5   |
| H A-1 LOCATION (X).....:          | 0.57      | 14.4  |
| J A-1 LOCATION (Y).....:          | 0.44      | 11.2  |
| K WELL TO WELL SPACING.....:      | 0.36      | 9.0   |

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



### CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

|      |              |         |     |              |      |
|------|--------------|---------|-----|--------------|------|
| WABF | MCC          | APPS    | PVC | PSS 120023XX | PR   |
| REV  | RD R 2/20/02 | 2/20/02 |     |              | 0202 |

ROBOTIC RESERVOIR, PP, CONVOLUTED BOTTOM