


CATALOG NUMBER	DESCRIPTION
1200-2300	CONVOLUTED BOTTOM, PP, NON-STERILE
1200-2301	CONVOLUTED BOTTOM, PP, STERILE

CATALOG NUMBER.....:	SEE ABOVE	
NOM. CAP (mL).....:	300	
MATERIAL.....:	PP	
	INCH	MM
A (BASE LENGTH).....:	5.03	127.8
B (BASE WIDTH).....:	3.365	85.5
C (TOP LENGTH).....:	4.88	123.8
D (TOP WIDTH).....:	3.22	81.7
E (STEP HEIGHT).....:	.10	2.5
F (OVERALL HEIGHT).....:	1.57	39.9
G (BOTTOM OF WELL DISTANCE).....:	.10	2.5
H A-1 LOCATION (X).....:	0.57	14.4
J A-1 LOCATION (Y).....:	0.44	11.2
K WELL TO WELL SPACING.....:	0.36	9.0

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



### CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

WABF	MCC	APPS	PVC	PSS 120023XX	PR
REV	RD R 2/20/02	2/20/02			0202

ROBOTIC RESERVOIR, PP, CONVOLUTED BOTTOM