




CATALOG NO.	5970-OX13	5970-OX16	5970-OX20
TEST TUBE SIZE (mm)	13	16	20
MATERIAL	RESMER	RESMER	RESMER
NOM. WEIGHT (gns)	120	200	190
	INCH	MM	INCH
A (LENGTH)	7.90±0.07	200.7±1.8	9.80±0.07
B (WIDTH)	4.06±0.05	103.1±1.3	4.98±0.05
C (HEIGHT)	2.30±0.05	58.4±1.3	2.84±0.05
D (SQUARE)	0.54±0.01	13.7±0.3	0.69±0.01
HOLE ARRAY	6 x 12	6 x 12	4 x 10

CATALOG NO.	5970-OX25	5970-OX30
TEST TUBE SIZE (mm)	25	30
MATERIAL	RESMER	RESMER
NOM. WEIGHT (gns)	250	220
	INCH	MM
A (LENGTH)	11.80±0.07	299.7±1.8
B (WIDTH)	4.84±0.05	122.9±1.3
C (HEIGHT)	3.75±0.05	95.3±1.3
D (SQUARE)	1.03±0.01	26.2±0.3
HOLE ARRAY	4 x 10	3 x 8

COLOR	CATALOG NUMBER			
WHITE	5970-0013	5970-0016	5970-0020	5970-0025
ORANGE	5970-0113	5970-0116	5970-0120	5970-0125
YELLOW	5970-0213	5970-0216	5970-0220	5970-0225
BLUE	5970-0313	5970-0316	5970-0320	5970-0325
GREEN	5970-0413	5970-0416	5970-0420	5970-0425
RED	5970-0513	5970-0516	5970-0520	5970-0525

- NOTES:
- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
 - THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Pancosmo Creek Drive
Rochester, NY 14625-2385, USA

MFG: NC
CR'G: GRP
APP'D: PVC
2/26/04
DESCRIPTION:

REV: 0204

PSS 00-5970-0XXX

UNWIRE TEST TUBE RACKS

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.