




CATALOG NO.	5970-OX13	5970-OX16	5970-OX20
TEST TUBE SIZE (mm)	13	16	20
MATERIAL	RESMER	RESMER	RESMER
NOM. WEIGHT (gns)	120	200	190
	INCH	MM	INCH
A (LENGTH)	7.90±0.07	200.7±1.8	9.80±0.07
B (WIDTH)	4.06±0.05	103.1±1.3	4.98±0.05
C (HEIGHT)	2.30±0.05	58.4±1.3	2.84±0.05
D (SQUARE)	0.54±0.01	13.7±0.3	0.69±0.01
HOLE ARRAY	6 x 12	6 x 12	4 x 10

CATALOG NO.	5970-OX25	5970-OX30
TEST TUBE SIZE (mm)	25	30
MATERIAL	RESMER	RESMER
NOM. WEIGHT (gns)	250	220
	INCH	MM
A (LENGTH)	11.80±0.07	299.7±1.8
B (WIDTH)	4.84±0.05	122.9±1.3
C (HEIGHT)	3.75±0.05	95.3±1.3
D (SQUARE)	1.03±0.01	26.2±0.3
HOLE ARRAY	4 x 10	3 x 8

COLOR	CATALOG NUMBER			
WHITE	5970-0013	5970-0016	5970-0020	5970-0025
ORANGE	5970-0113	5970-0116	5970-0120	5970-0125
YELLOW	5970-0213	5970-0216	5970-0220	5970-0225
BLUE	5970-0313	5970-0316	5970-0320	5970-0325
GREEN	5970-0413	5970-0416	5970-0420	5970-0425
RED	5970-0513	5970-0516	5970-0520	5970-0525

- NOTES:
- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
  - THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Pancosmo Creek Drive  
Rochester, NY 14625-2385, USA

MFG	NC	APP'D	PVC
CHK'D	GRP	2/26/04	
DESCRIPTION:			

REV: 0204

PSS 00-5970-0XXX

UNWIRE TEST TUBE RACKS