


CATALOG NO.	NOMINAL TUBE SIZE, (mm)				MATERIAL	
	13	16	20	25		30
005972						
SEE BELOW						
SIZE CODE	7.5		15		18	
MAX. TUBE CAP. (nL)	62		118		108	
NOM. WEIGHT (gms)						
	INCHES	MM	INCHES	MM	INCHES	MM
A (LENGTH)	4.06 ±0.05	103.1 ±1.3	4.99 ±0.06	126.7 ±1.5	5.06 ±0.06	128.5 ±1.5
B (WIDTH)	4.04 ±0.05	102.6 ±1.3	4.97 ±0.06	126.2 ±1.5	4.06 ±0.05	103.1 ±1.3
C (HEIGHT)	2.31 ±0.04	58.7 ±1.0	2.80 ±0.04	71.1 ±1.0	3.34 ±0.04	84.8 ±1.0
D (SQUARE)	0.56 ±0.01	14.2 ±0.3	0.68 ±0.02	17.3 ±0.5	0.85 ±0.02	21.6 ±0.5
HOLE ARRAY	6 x 6		6 x 6		4 x 5	
COLOR	SIZE CODE					
WHITE	-0013		-0016		-0020	
BLUE	-0313		-0316		-0320	
GREEN	-0413		-0416		-0420	
RED	-0513		-0516		-0520	
SIZE CODE	SEE BELOW					
MAX. TUBE CAP. (nL)	25		30		108	
NOM. WEIGHT (gms)	120		108		108	
	INCHES	MM	INCHES	MM	INCHES	MM
A (LENGTH)	4.83 ±0.06	122.7 ±1.5	4.30 ±0.05	109.2 ±1.3	4.83 ±0.06	122.7 ±1.5
B (WIDTH)	4.82 ±0.06	122.4 ±1.5	4.30 ±0.05	109.2 ±1.3	4.82 ±0.06	122.4 ±1.5
C (HEIGHT)	2.95 ±0.04	74.9 ±1.0	3.36 ±0.04	85.3 ±1.0	2.95 ±0.04	74.9 ±1.0
D (SQUARE)	1.03 ±0.02	26.2 ±0.5	1.25 ±0.02	31.8 ±0.5	1.03 ±0.02	26.2 ±0.5
HOLE ARRAY	4 x 4		3 x 3		3 x 3	
COLOR	SIZE CODE					
WHITE	-0025		-0030		-0030	
BLUE	-0325		-0330		-0330	
GREEN	-		-0430		-0430	
RED	-		-0530		-0530	

- NOTES:
1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
 2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



CONFIDENTIAL
 THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
 75 Pancosmo Creek Drive
 Rochester, NY 14625-2385, USA

NAME	MC	APP'D PVC	PSS 00-5972-0XXX-1	REV.	0204
DATE	PVC	2/26/04			
DESCRIPTION:					

LINWIRE HALF-RACKS

SPECS ON BACK

NOT TO SCALE