


CATALOG NO.	NOMINAL CAPACITY, (oz)							MATERIALS	
								BOTTLE	CLOSURE
002016	1	2	4	8	16	32	-	PP	PP
002018	1	2	4	8	16	32	-	HOPE	PP

SIZE CODE.....	-0030		-0060		-0125	
SIZE (NOM VOL mL/oz)	30/1		60/2		125/4	
NOM. BRIM CAP. (mL):	45		80		180	
NECK FINISH (mm)....	20-415		24-415		38-430	
THREAD (T.P.I.).....	8		8		6	
NOM. WEIGHT W/CLOSURE (gms).....	10		14		40	
	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.).....	0.55 ±0.02	14.0 ±0.5	0.70 ±0.02	17.8 ±0.5	1.100 ±0.025	27.9 ±0.64
B (HT. W/CLOSURE)...	2.50 ±0.07	63.5 ±1.8	3.30 ±0.08	83.8 ±2.0	4.32 ±0.08	109.7 ±2.0
C (HT. W/O CLOSURE):	2.43 ±0.07	61.7 ±1.8	3.21 ±0.08	81.5 ±2.0	4.15 ±0.08	105.4 ±2.0
D (SQUARE).....	1.52 ±0.06	38.6 ±1.5	1.61 ±0.06	40.9 ±1.5	2.12 ±0.07	53.8 ±1.8
W.T.(NOM. WALL THK):	0.04	1.0	0.04	1.0	0.06	1.5

SIZE CODE.....	-0250		-0500		-1000	
SIZE (NOM VOL mL/oz)	250/8		500/16		1000/32	
NOM. BRIM CAP. (mL):	330		620		1230	
NECK FINISH (mm)....	38-430		38-430		38-430	
THREAD (T.P.I.).....	6		6		6	
NOM. WEIGHT W/CLOSURE (gms).....	50		75		120	
	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.).....	1.100 ±0.025	27.9 ±0.64	1.100 ±0.025	27.9 ±0.64	1.100 ±0.025	27.9 ±0.64
B (HT. W/CLOSURE)...	5.74 ±0.08	145.8 ±2.0	6.98 ±0.09	177.3 ±2.3	8.65 ±0.10	219.7 ±2.5
C (HT. W/O CLOSURE):	5.58 ±0.08	141.7 ±2.0	6.81 ±0.09	173.0 ±2.3	8.46 ±0.10	214.9 ±2.5
D (SQUARE).....	2.41 ±0.07	61.2 ±1.8	2.93 ±0.08	74.4 ±2.0	3.70 ±0.08	94.0 ±2.0
W.T.(NOM. WALL THK):	0.06	1.5	0.065	1.7	0.05	1.3

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE CK'D	MC <sup>2</sup>	APP'D	PSS 00-201X-XXXX-1	REV. 1103
--------------	-----------------	-------	--------------------	--------------

DESCRIPTION:  
SQUARE, NARROW-MOUTH  
(GRADUATED) LAB QUALITY BOTTLES