


CATALOG NO.	NOMINAL CAPACITY, (oz)							MATERIALS	
								BOTTLE	CLOSURE
002002	1	2	4	6	8	16	32	HOPE	PP
002003	1	2	4	—	8	16	32	LOPE	PP
002004	1	2	4	—	8	16	32	AMBER HOPE	AMBER PP
002006	1	2	4	—	8	16	32	PP	PP

SIZE CODE.....	-0001		-0002		-0004		-0006	
SIZE (NOM VOL mL/oz)	30/1		60/2		125/4		175/6	
NOM. BRIM CAP. (mL):	34		64		140		185	
NECK FINISH (mm)....	20-415		20-415		24-415		24-415	
THREAD (T.P.I.).....	8		8		8		8	
NOM. WEIGHT W/CLOSURE (gms).....	8		12		20		25	
	INCH	MM	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.).....	0.54 ±0.02	13.7 ±0.5	0.54 ±0.02	13.7 ±0.5	0.69 ±0.02	17.5 ±0.5	0.69 ±0.02	17.5 ±0.5
B (HT. W/CLOSURE)....	2.40 ±0.07	61.0 ±1.8	3.36 ±0.08	85.3 ±2.0	3.97 ±0.08	100.8 ±2.0	4.90 ±0.08	124.5 ±2.0
C (HT. W/O CLOSURE):	2.30 ±0.07	58.4 ±1.8	3.25 ±0.08	82.6 ±2.0	3.88 ±0.08	98.6 ±2.0	4.81 ±0.08	122.2 ±2.0
D (DIA).....	1.35 ±0.05	34.3 ±1.3	1.52 ±0.06	38.6 ±1.5	1.98 ±0.06	50.3 ±1.5	2.13 ±0.06	54.1 ±1.5
W.T.(NOM. WALL THK):	0.035	0.90	0.035	0.90	0.040	1.02	0.040	1.02

SIZE CODE.....	-0008		-0016		-9016		-0032	
SIZE (NOM VOL mL/oz)	250/8		500/16		500/16		1000/32	
NOM. BRIM CAP. (mL):	285		525		525		1045	
NECK FINISH (mm)....	24-415		28-415		38-430		38-430	
THREAD (T.P.I.).....	8		6		6		6	
NOM. WEIGHT W/CLOSURE (gms).....	32		58		52		100	
	INCH	MM	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.).....	0.69 ±0.02	17.5 ±0.5	0.84 ±0.02	21.3 ±0.5	1.08 ±0.03	27.4 ±0.8	1.08 ±0.03	27.4 ±0.8
B (HT. W/CLOSURE)....	5.22 ±0.08	132.6 ±2.0	6.71 ±0.09	170.4 ±2.3	6.71 ±0.09	170.4 ±2.3	8.50 ±0.10	215.9 ±2.5
C (HT. W/O CLOSURE):	5.13 ±0.08	130.3 ±2.0	6.59 ±0.09	167.4 ±2.3	6.59 ±0.09	167.4 ±2.3	8.35 ±0.10	212.1 ±2.5
D (DIA).....	2.40 ±0.07	61.0 ±1.8	2.86 ±0.08	72.6 ±2.0	2.86 ±0.08	72.6 ±2.0	3.60 ±0.08	91.4 ±2.0
W.T.(NOM. WALL THK):	0.050	1.27	0.050	1.27	0.050	1.27	0.050	1.27

- NOTES:
- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
 - THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher SCIENTIFIC		Thermo Fisher Scientific 75 Panorama Creek Drive Rochester, NY 14625-2385, USA	
MADE JJ/MG PVC	APP'D PVC 12/12/00	PSS 00-200X-XXXX-2	REV. 1200
DESCRIPTION: ROUND NARROW MOUTH LAB QUALITY BOTTLES			