

# Thermo Scientific Nalgene Economy Bottles

## Narrow Mouth Bottles



It's easy to choose high-quality lab bottles and stay within budget. Nalgene® Economy Bottles are lightweight versions of the original Nalgene Lab Quality Bottles.

These bottles are suitable for light-duty, general purpose laboratory applications such as containing samples and daily working solutions like buffers and bench reagents. Like most Nalgene bottles, the Economy Bottles are guaranteed leakproof and are made in our ISO-certified USA manufacturing facility.

## Wide Mouth Bottles



Available in chemical-resistant high density polyethylene (HDPE) or autoclavable polypropylene copolymer (PPCO), they are selected for a wide variety of laboratory applications. Priced below low-quality plastic or glass bottle brands, Nalgene Economy Bottles are the economical, quality choice for lighter-duty lab use.

## Anatomy of a Nalgene Economy Bottle

### Closure

is one-piece linerless.

The patented closure and bottle system, work together to supply a guaranteed leakproof system.

### Shrink Ring

at the neck of most NALGENE bottles is molded so the inside surface is smooth, minimising content entrapment.

### Uniform Walls

are resistant to splitting or puncturing.

Seal Ring is molded inside the closure. It fits tightly against the beveled inner edge (chamfer) of the bottle neck. No closure liner to wear, crease or cause contamination.

Threads on bottles and closures have continuous, straight-shouldered semi-buttress threads, not low-quality round threads.

Bottom has curved inner corners for easy cleaning. Stable base has permanent, molded-in resin code and volume.



The patented Nalgene bottle & closure system – the only *guaranteed* leakproof system available.