


CATALOG NO.	NOMINAL CAPACITY, (mL/oz)		MATERIALS
001223	1000/32	3000/96	PMP

SIZE CODE.....:	-1000		-3000	
SIZE (NOM VOL mL/oz):	1000/32		3000/96	
NOM. WEIGHT (gms)...	190		330	
	IN	MM	IN	MM
A	6.5	164	8.9	226
B (DIA).....	3.8	97	6.3	160
H	7.2	183	7.6	193

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Ponorona Creek Drive
Rochester, NY 14625-2385, USA

MADE	NS	APP'D	PVC	PSS 00-1223-XXXX-1	REV.	—
CHK'D	PVC		2/7/91			

DESCRIPTION:

GRADUATED BEAKERS WITH HANDLE