


CATALOG NO.	NOMINAL CAPACITY, (oz)							MATERIALS	
								BOTTLE	CLOSURE
001600	1	2	4	8	16	32	64	TEFLON FEP	TEFZEL ETFE
001620	-	-	4	8	16	32	-	BLACK TEFLON FEP	BLACK TEFZEL ETFE

SIZE CODE.....	-0001	-0002	-0004	-0008				
SIZE (NOM VOL mL/oz)	30/1	60/2	125/4	250/8				
NOM. BRIM CAP. (mL)	35	70	140	280				
NECK FINISH (mm)...	20-415	20-415	24-415	24-415				
THREAD (T.P.I.).....	8	8	8	8				
NOM. WEIGHT W/CLOSURE (gms).....	20	22	40	60				
	INCH	MM	INCH	MM	INCH	MM	INCH	MM
A (BTLE. NECK I.D.)	0.55 ±0.03	14.0 ±0.8	0.55 ±0.03	14.0 ±0.8	0.68 ±0.03	17.3 ±0.8	0.68 ±0.03	17.3 ±0.8
B (HT. W/CLOSURE)...	2.95 ±0.06	74.9 ±1.5	3.32 ±0.07	84.3 ±1.8	4.52 ±0.08	114.8 ±2.0	5.28 ±0.08	134.1 ±2.0
C (HT. W/O CLOSURE)...	2.86 ±0.06	72.6 ±1.5	3.23 ±0.07	82.0 ±1.8	4.40 ±0.08	111.8 ±2.0	5.16 ±0.08	131.1 ±2.0
D (DIA).....	1.25 ±0.03	31.8 ±0.8	1.55 ±0.05	39.4 ±1.3	1.80 ±0.05	45.7 ±1.3	2.36 ±0.06	59.9 ±1.5
W.T.(NOM. WALL THK):	0.040	1.02	0.040	1.02	0.040	1.02	0.030	0.76

SIZE CODE.....	-0016	-0032	-0064			
SIZE (NOM VOL mL/oz)	500/16	1000/32	2000/64			
NOM. BRIM CAP. (mL)	520	1040	2100			
NECK FINISH (mm)...	28-415	38-415	38-430			
THREAD (T.P.I.).....	6	6	6			
NOM. WEIGHT W/CLOSURE (gms).....	100	160	340			
	INCH	MM	INCH	MM	INCH	MM
A (BTLE. NECK I.D.)	0.80 ±0.03	20.3 ±0.8	1.04 ±0.03	26.4 ±0.8	0.96 ±0.04	24.4 ±1.0
B (HT. W/CLOSURE)...	6.53 ±0.09	165.9 ±2.3	7.98 ±0.09	202.7 ±2.3	9.66 ±0.10	245.4 ±2.5
C (HT. W/O CLOSURE)...	6.38 ±0.09	162.1 ±2.3	7.78 ±0.09	197.6 ±2.3	9.42 ±0.10	239.3 ±2.5
D (DIA).....	2.86 ±0.06	72.6 ±1.5	3.55 ±0.07	90.2 ±1.8	4.75 ±0.08	120.7 ±2.0
W.T.(NOM. WALL THK):	0.040	1.02	0.025	0.64	0.050	1.27

NOTES:

- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
- THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	JGJ	APP'D	PVC	PSS 00-16XX-XXXX-1	REV.	0898
CHK'D	PVC	DATE	8/26/98			
DESCRIPTION:				NARROW-MOUTH ROUND LAB QUALITY BOTTLES		