


CATALOG NO.	NOMINAL CAPACITY, (oz)						MATERIALS	
							BOTTLE	CLOSURE
001630	1	2	4	8	16	32	TEFLON PFA	TEFLON PFA

SIZE CODE.....	-0001		-0002		-0004	
SIZE (NOM VOL mL/oz)	30/1		60/2		125/4	
NOM. BRIM CAP. (mL)	35		70		150	
NECK FINISH (mm)....	20-415		20-415		38-430	
THREAD (T.P.I.).....	8		8		6	
NOM. WEIGHT W/CLOSURE (gns).....	20		25		60	
	INCH	MM	INCH	MM	INCH	MM
A (BTLE. NECK I.D.):	0.55 ± 0.03	13.9±0.8	0.55 ± 0.03	14.0±0.8	0.96 ± 0.04	24.4±1.0
B (HT. W/CLOSURE)...	3.00 ± 0.07	76.2±1.8	3.40 ± 0.08	86.4±2.0	5.00 ± 0.09	127.0±2.3
C (HT. W/O CLOSURE):	2.85 ± 0.07	72.4±1.8	3.25 ± 0.08	82.6±2.0	4.80 ± 0.09	121.9±2.3
D (DIA).....	1.25 ± 0.04	31.8±1.0	1.55 ± 0.05	39.4±1.3	1.80 ± 0.06	45.7±1.5
W.T.(NOM. WALL THK):	0.040	1.0	0.040	1.0	0.040	1.0

SIZE CODE.....	-0008		-0016		-0032	
SIZE (NOM VOL mL/oz)	250/8		500/16		1000/32	
NOM. BRIM CAP. (mL)	280		520		1040	
NECK FINISH (mm)....	38-430		38-430		38-430	
THREAD (T.P.I.).....	6		6		6	
NOM. WEIGHT W/CLOSURE (gns).....	85		125		180	
	INCH	MM	INCH	MM	INCH	MM
A (BTLE. NECK I.D.):	0.96 ± 0.04	24.4±1.0	0.96 ± 0.04	24.4±1.0	0.96 ± 0.04	24.4±1.0
B (HT. W/CLOSURE)...	5.75 ± 0.09	146.1±2.3	7.10 ± 0.10	180.3±2.5	8.50 ± 0.10	215.9±2.5
C (HT. W/O CLOSURE):	5.50 ± 0.09	139.7±2.3	6.90 ± 0.10	175.3±2.5	8.25 ± 0.10	209.6±2.5
D (DIA).....	2.30 ± 0.07	58.4±1.8	2.80 ± 0.07	71.1±1.8	3.55 ± 0.08	90.2±2.0
W.T.(NOM. WALL THK):	0.040	1.0	0.040	1.0	0.040	1.0

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MCC	APP'D	GRP	PSS 00-1630-XXXX-1	REV.	0907
DATE	9/24/07	DATE	9/24/07			

DESCRIPTION:

NARROW-MOUTH ROUND  
LAB QUALITY BOTTLES