


CATALOG NUMBER.....	14100-0002		14100-0010		14100-0015		14100-0020	
	2		10		15		20	
NOM. CAP (GAL).....	LLDPE		LLDPE		LLDPE		LLDPE	
MATERIAL.....	INCH	MM	INCH	MM	INCH	MM	INCH	MM
A (LENGTH).....	10 1/2	267	14 9/16	370	20 11/16	525	20 11/16	525
B (WIDTH).....	10 1/2	267	14 9/16	370	14 1/2	368	14 5/8	371
C (HT. W/O COVER)....	9 1/4	235	14 5/16	364	13 9/16	235	19 13/16	503
D (LENGTH).....	8 3/16	208	12 1/4	311	18 7/16	344	18 7/16	468
E (WIDTH).....	8 3/16	208	12 1/4	311	12 1/2	318	12 1/2	318
F (LENGTH COVER)....	11 3/16	284	15 1/16	383	21 5/16	541	21 5/16	541
G (WIDTH COVER)....	11 3/16	284	15 1/16	383	15 1/4	387	15 1/4	387
H (HT. COVER).....	1 1/4	32	1 1/4	32	1 1/4	32	1 1/4	32
W.T. (NOM. WALL THK.)	5/32	4	5/32	4	5/32	4	5/32	4

CATALOG NUMBER.....	14100-0040		14100-0045		14100-0065	
	40		45		65	
NOM. CAP (GAL).....	LLDPE		LLDPE		LLDPE	
MATERIAL.....	INCH	MM	INCH	MM	INCH	MM
A (LENGTH).....	26 5/8	676	26 5/8	676	26 11/16	678
B (WIDTH).....	14 1/2	368	20 11/16	525	20 9/16	522
C (HT. W/O COVER)....	13 1/4	337	19 9/16	483	25 3/16	640
D (LENGTH).....	24 3/8	619	24 3/16	614	24 7/16	621
E (WIDTH).....	12 1/2	318	18 7/16	468	18 5/16	465
F (LENGTH COVER)....	27 1/8	689	27 7/16	697	27 7/16	697
G (WIDTH COVER)....	15 1/8	384	21 1/8	537	21 1/8	537
H (HT. COVER).....	1 1/4	32	1 1/4	32	1 1/4	32
W.T. (NOM. WALL THK.)	5/32	4	5/32	4	5/32	4

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

DATE	PB	DATE	GRP
BY	RCB	8/20/13	

PSS 14100-00XX

REV. 0813

TANKS WITH COVERS