


CATALOG NUMBER.....	2135-3803	2135-5303	2135-8303
CLOSURE.....	PP		
MATERIAL INSERT.....	SILICONE - PLATINUM CURED		
TUBING.....	SILICONE - PLATINUM CURED		
CLOSURE SIZE (mm).....	38-430	53B	83B
TUBING SIZES (COMBINATION).....	(2) 1/8 I.D. (1) 1/4 I.D.	(2) 1/4 I.D. (1) 1/8 I.D.	(2) 1/4 I.D. (1) 3/8 I.D.
A DIA.....	INCH	MM	INCH
	1.66±.03	42.2±0.8	2.58±.02
B (CLOSURE HT).....	1.17±.03	29.7±0.8	1.65±.02
H (INSERT HT).....	0.42±.03	10.7±0.8	0.67±.03
TUBING LENGTHS (ABOVE INSERT).....	12	305	12
	1/8 I.D.	406	6
TUBING LENGTHS (BELOW INSERT)	6	152	20
	3/8 I.D.	—	—
NOM WT. (gms).....	80	140	345



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL
 THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
 75 Pennington Creek Drive
 Rochester, NY 14625-2385, USA

REV	0501
DATE	5/15/01
DESCRIPTION	FLEXIBLE TOP WORKS™ SYSTEM SCREW CLOSURES WITH SILICONE INSERTS AND 3 PORTS
PART NO.	PSS 00-2135-XX03

NOTES:
 1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
 2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.