


CATALOG NO.	NOMINAL CAPACITY, (oz)					MATERIALS	
						BOTTLE	CLOSURE
002097	8	16	32	64	1 GAL	FLPE	FLPP

SIZE CODE.....:	-0008		-0016		-0032	
SIZE (NOM VOL mL/oz):	250/8		500/16		1000/32	
NOM. BRIM CAP. (mL):	285		525		1045	
NECK FINISH (mm)....:	24-415		28-415		38-430	
THREAD (T.P.I.).....:	8		6		6	
NOM. WEIGHT W/CLOSURE (gms).....:	32		58		100	
	INCH	MM	INCH	MM	INCH	MM
A (BTL NECK I.D.)...:	0.69 ±0.02	17.5 ±0.5	0.84 ±0.02	21.3± 0.5	1.08 ±0.03	27.4 ±0.8
B (HT. W/CLOSURE)...:	5.22 ±0.08	132.6 ±2.0	6.71 ±0.09	170.4 ±2.3	8.50 ±0.10	215.9 ±2.5
C (HT. W/D CLOSURE):	5.13 ±0.08	130.3 ±2.0	6.59 ±0.09	167.4 ±2.3	8.35 ±0.10	212.1 ±2.5
D (DIA).....:	2.40 ±0.07	61.0 ±1.8	2.86 ±0.08	72.6 ±2.0	3.60 ±0.08	91.4 ±2.0
W.T.(NOM. WALL THK):	0.050	1.27	0.050	1.27	0.050	1.27

SIZE CODE.....:	-0005		-0010	
SIZE (NOM VOL mL/oz):	2000/64		4000/1gal	
NOM. BRIM CAP. (mL):	2100		4160	
NECK FINISH (mm)....:	38-430		38-430	
THREAD (T.P.I.).....:	6		6	
NOM. WEIGHT W/CLOSURE (gms).....:	160		370	
	INCH	MM	INCH	MM
A (BTL NECK I.D.)...:	1.08 ±0.03	27.4 ±0.8	1.02 ±0.03	25.9 ±0.8
B (HT. W/CLOSURE)...:	9.85 ±0.12	250.2 ±3.0	13.14 ±0.16	333.8 ±4.1
C (HT. W/D CLOSURE):	9.68 ±0.12	245.9 ±3.0	12.96 ±0.16	329.2 ±4.1
D (DIA).....:	4.74 ±0.08	120.4 ±2.0	6.02 ±0.10	152.9 ±2.5
W.T.(NOM. WALL THK):	0.050	1.27	0.050	1.27

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE IN	MG	APP'D	PVC	PSS 00-2097-XXXX-1	REV.
CK'D	PVC		4/28/99		0499

DESCRIPTION:  
FLUORINATED NARROW-MOUTH BOTTLES