


CATALOG NO.	NOMINAL CAPACITY ( GAL )		MATERIALS			
			CARBOY	CLOSURE	GASKET IN CLOSURE	SPIGOT (WITH TEFLON O-RING)
002317	2 1/2	5	TRANSPARENT PC	WHITE PP	TPE (KRATON)	PP

SIZE CODE.....:	-0020		-0050	
SIZE (NOM VOL L/GAL):	10 / 2.5		20 / 5.0	
NOM. BRIM VOL (L)....:	12.5		24	
NECK FINISH (mm)....:	83B		83B	
THREAD (T.P.I.).....:	2		2	
NOM. WEIGHT W/ CLOSURE(gms).....:	900		1537	
	INCHES	MM	INCHES	MM
A (NECK I.D.).....:	2.70 <sup>+0.12</sup> <sub>-0.16</sub>	68.6 <sup>+3.0</sup> <sub>-4.1</sub>	2.70 <sup>+0.12</sup> <sub>-0.16</sub>	68.6 <sup>+3.0</sup> <sub>-4.1</sub>
B (HT. W/ CLOSURE)...:	15.5 ± 0.3	394 ± 8	21.1 ± 0.4	536 ± 10
C (HT. W/O CLOSURE)..:	15.0 ± 0.3	381 ± 8	20.6 ± 0.4	523 ± 10
D (DIA).....:	9.96 ± 0.16	253.0 ± 4.1	11.3 ± 0.2	287 ± 5
W.T.(NOM WALL THK)...:	0.07	1.8	0.07	1.8

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher  
SCIENTIFIC**

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MG	APP'D	PSS 00-2317-XXXX-1	REV.
CK'D	KR	GRP		0307

DESCRIPTION:  
ROUND (CLEARBOY™) CARBOY