







CAT. NO.....:	455-0150	455-0250	455-0500	455-1000				
NOM. VOLUME (mL)...:	150	250	500	1000				
MATERIALS	RECIEVER...:	PS						
	STORAGE CLOSURE...:	WHITE HDPE						
WT W/ CLOSURE (gm)...:	50		62		95		145	
	INCHES	MM	INCHES	MM	INCHES	MM	INCHES	MM
A DIA (NECK I.D.)...:	1.47	37.3	1.47	37.3	1.47	37.3	1.47	37.3
B (HT W/O CLOSURE)...:	3.33	84.6	4.25	108.0	5.50	139.7	7.14	181.4
C (HT W/ CLOSURE)...:	3.50	88.9	4.42	112.3	5.67	144.0	7.32	185.9
D DIA (O.D.)...:	3.50	88.9	3.50	88.9	3.86	98.0	4.48	113.8
E DIA (CLOSURE)...:	1.96	49.8	1.96	49.8	1.96	49.8	1.96	49.8
F (CLOSURE HT)...:	1.12	28.4	1.12	28.4	1.12	28.4	1.12	28.4



NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MG/MC	APP'D	PVC	PSS 00-0455-XXXX-1	REV.	0699
CK'D	PVC		6/25/99			

DESCRIPTION:

FILTER UNIT RECEIVERS