

Nalgene Centrifuge Bottles with Sealing Closure

polypropylene copolymer, polypropylene closure

Thermo Scientific™ Nalgene™ centrifuge bottles are ideal for large-volume cell harvesting, pelleting and protein purification.

Molded of transparent PC, these translucent bottles are autoclavable and have excellent mechanical strength properties.

details

- For low- to moderate-speed centrifugation of biological and chemical samples
- Molded clear PC for easy visualization of samples
- PP sealing cap closure with silicone gasket, for added leakproof¹ assurance at maximum speeds or when spinning hazardous materials

Alert: When spinning 250 mL bottle in a Sorvall™ GSA rotor, use of the Nalgene centrifuge bottle adapter (DS3125-0250) is required.

Note: For proper performance bottles must be filled to 80% of total capacity



Nalgene Centrifuge Bottles with Sealing Closures

Cat. No.	Capacity, mL	OD x H, mm	Max Force, xg	No. per Pack/Case
3140-0250	250	61,8 x 134,4	27,500	4/36
3140-0500	450	69,5 x 159,5	13,700	4/24