

Nalgene Centrifuge Bottles with Sealing Closure

polypropylene copolymer, polypropylene closure

Thermo Scientific™ Nalgene™ centrifuge bottles are ideal for large-volume cell harvesting, pelleting and protein purification.

Molded of rigid and robust PPCO, these translucent bottles are autoclavable and have excellent chemical resistance properties.

details

- For low- to moderate-speed centrifugation of biological and chemical samples
- Bottle includes PP sealing cap assembly with silicone gasket for added leakproof assurance at maximum speeds or when spinning hazardous materials

Requires: When spinning 250 mL bottle in a Sorvall GSA rotor, use of the Nalgene centrifuge bottle adapter (DS3125-0250) is required.

Note: For proper performance bottles must be filled to 80% of total capacity

Nalgene Centrifuge Bottles with Sealing Closure

Cat. No.	Capacity, mL	OD x H, mm	Max Force, xg	No. per Pack/Case
3141-0250	250	61.1 x 132.8	27,500	4/36
3141-0500	450	69.5 x 159.3	13,700	4/24