

| CATALOG NO. | NOMINAL CAPACITY (GAL.) | | | MATERIALS | | |
|-------------|-------------------------|---|---|-----------------|--------------------|-----------|
| | | | | BOTTLE | CLOSURE | SEAL RING |
| 002202 | 1/2 | 1 | 2 | LOPE | PP | — |
| 002203 | 1/2 | 1 | 2 | PP | PP | TPE*** |
| 002204 | 1/2 | 1 | 2 | OPAQUE AMBER PP | OPAQUE AMBER PP ** | TPE*** |
| DS2205 | 1/2 | 1 | 2 | PC | PP | TPE*** |


NOTES:

1. SIZE CODES -0210 AND -0250 ARE FOR CATALOG NUMBER DS2205.
- * 2. 1 AND 2 GALLON SIZES HAVE MOLDED-IN SHOULDER LOOP (FOR ATTACHING IDENTIFICATION TAG) AS SHOWN.
- ** 3. THE CLOSURE ON CATALOG NUMBER 002204-0020 IS WHITE PP.
- *** 4. SEAL RING APPLIES ONLY TO -0020 SIZE FOR THESE CATALOGS.

| SIZE CODE..... | -0005/-0210 | | -0250 | | -0010* | | -0020* | |
|---------------------------------|-------------|------------|------------|------------|------------|------------|-----------|------------|
| SIZE (NOM VOL L/GAL) | 2/0.5 | | 2.5/0.5 | | 4/1 | | 8/2 | |
| NOM. BRIM CAP. (mL): | 2100 | | 2275 | | 4300 | | 8000 | |
| NECK FINISH (mm).... | 38-430 | | 38-430 | | 38-430 | | 53B | |
| THREAD (T.P.I.).... | 6 | | 6 | | 6 | | 4 | |
| NOM. WEIGHT W/CLOSURE (gms).... | 180 | | 195 | | 370 | | 700 | |
| | INCHES | MM | INCHES | MM | INCHES | MM | INCHES | MM |
| A (NECK I.D.)..... | 1.00 ±.05 | 25.4 ±1.3 | 1.00 ±.05 | 25.4 ±1.3 | 1.00 ±.05 | 25.4 ±1.3 | 1.46 ±.10 | 37.1 ±2.5 |
| B (HT. W/CLOSURE).. | 9.80 ±.12 | 248.9 ±3.0 | 11.60 ±.14 | 294.6 ±3.6 | 13.10 ±.15 | 332.7 ±3.8 | 16.1 ±.3 | 409 ±8 |
| C (HT. W/O CLOSURE): | 9.60 ±.12 | 243.8 ±3.0 | 11.40 ±.14 | 289.6 ±3.6 | 12.90 ±.15 | 327.7 ±3.8 | 15.9 ±.2 | 404 ±5 |
| D (DIA)..... | 4.75 ±.10 | 120.7 ±2.5 | 4.75 ±.10 | 120.7 ±2.5 | 6.00 ±.12 | 152.4 ±3.0 | 7.60 ±.13 | 193.0 ±3.3 |
| W.T.(NOM WALL THK).. | 0.07 | 1.8 | 0.08 | 2.0 | 0.10 | 2.5 | 0.12 | 3.0 |

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

| | | | |
|--------------|-----------|-----------------------------------|-----------|
| MADE MC | APP'D PVC | PSS 00-220X-XXXX-1 | REV. 0403 |
| CHK'D PVC | 4/24/03 | | |
| DESCRIPTION: | | LARGE, NARROW-MOUTH ROUND BOTTLES | |