

CATALOG NO.	NOMINAL CAPACITY (GAL.)			MATERIALS		
				BOTTLE	CLOSURE	SEAL RING
002202	1/2	1	2	LOPE	PP	—
002203	1/2	1	2	PP	PP	TPE***
002204	1/2	1	2	OPAQUE AMBER PP	OPAQUE AMBER PP **	TPE***
DS2205	1/2	1	2	PC	PP	TPE***


**NOTES:**

1. SIZE CODES -0210 AND -0250 ARE FOR CATALOG NUMBER DS2205.
- \* 2. 1 AND 2 GALLON SIZES HAVE MOLDED-IN SHOULDER LOOP (FOR ATTACHING IDENTIFICATION TAG) AS SHOWN.
- \*\* 3. THE CLOSURE ON CATALOG NUMBER 002204-0020 IS WHITE PP.
- \*\*\* 4. SEAL RING APPLIES ONLY TO -0020 SIZE FOR THESE CATALOGS.

SIZE CODE.....:	-0005/-0210	-0250	-0010*	-0020*
SIZE (NOM VOL L/GAL):	2/0.5	2.5/0.5	4/1	8/2
NOM. BRIM CAP. (mL):	2100	2275	4300	8000
NECK FINISH (mm):	38-430	38-430	38-430	53B
THREAD (T.P.I.):	6	6	6	4
NOM. WEIGHT W/CLOSURE (gms):	180	195	370	700
	INCHES	MM	INCHES	MM
A (NECK I.D.):	1.00 ± 0.05	25.4 ± 1.3	1.00 ± 0.05	25.4 ± 1.3
B (HT. W/CLOSURE):	9.80 ± 0.12	248.9 ± 3.0	11.60 ± 0.14	294.6 ± 3.6
C (HT. W/O CLOSURE):	9.60 ± 0.12	243.8 ± 3.0	11.40 ± 0.14	289.6 ± 3.6
D (DIA):	4.75 ± 0.10	120.7 ± 2.5	4.75 ± 0.10	120.7 ± 2.5
W.T. (NOM WALL THK):	0.07	1.8	0.08	2.0

**NOTES:**

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



**CONFIDENTIAL**

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher SCIENTIFIC**

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE MC	APP'D PVC	PSS 00-220X-XXXX-1	REV. 0403
CK'D PVC	4/24/03		
DESCRIPTION:		LARGE, NARROW-MOUTH ROUND BOTTLES	