

| CATALOG NO. | NOMINAL CAPACITY (GAL.) |   |   | MATERIALS       |                    |           |
|-------------|-------------------------|---|---|-----------------|--------------------|-----------|
|             |                         |   |   | BOTTLE          | CLOSURE            | SEAL RING |
| 002202      | 1/2                     | 1 | 2 | LDPE            | PP                 | —         |
| 002203      | 1/2                     | 1 | 2 | PP              | PP                 | TPE***    |
| 002204      | 1/2                     | 1 | 2 | OPAQUE AMBER PP | OPAQUE AMBER PP ** | TPE***    |
| DS2205      | 1/2                     | 1 | 2 | PC              | PP                 | TPE***    |


**NOTES:**

1. SIZE CODES -0210 AND -0250 ARE FOR CATALOG NUMBER DS2205.
- \* 2. 1 AND 2 GALLON SIZES HAVE MOLDED-IN SHOULDER LOOP (FOR ATTACHING IDENTIFICATION TAG) AS SHOWN.
- \*\* 3. THE CLOSURE ON CATALOG NUMBER 002204-0020 IS WHITE PP.
- \*\*\* 4. SEAL RING APPLIES ONLY TO -0020 SIZE FOR THESE CATALOGS.

| SIZE CODE.....                   | -0005/-0210 |            | -0250      |            | -0010*     |            | -0020*    |            |
|----------------------------------|-------------|------------|------------|------------|------------|------------|-----------|------------|
| SIZE (NOM VOL L/GAL)             | 2/0.5       |            | 2.5/0.5    |            | 4/1        |            | 8/2       |            |
| NOM. BRIM CAP. (mL):             | 2100        |            | 2275       |            | 4300       |            | 8000      |            |
| NECK FINISH (mm)....             | 38-430      |            | 38-430     |            | 38-430     |            | 53B       |            |
| THREAD (T.P.I.).....             | 6           |            | 6          |            | 6          |            | 4         |            |
| NOM. WEIGHT W/CLOSURE (gms)..... | 180         |            | 195        |            | 370        |            | 700       |            |
|                                  | INCHES      | MM         | INCHES     | MM         | INCHES     | MM         | INCHES    | MM         |
| A (NECK I.D.).....               | 1.00 ±.05   | 25.4 ±1.3  | 1.00 ±.05  | 25.4 ±1.3  | 1.00 ±.05  | 25.4 ±1.3  | 1.46 ±.10 | 37.1 ±2.5  |
| B (HT. W/CLOSURE)...             | 9.80 ±.12   | 248.9 ±3.0 | 11.60 ±.14 | 294.6 ±3.6 | 13.10 ±.15 | 332.7 ±3.8 | 16.1 ±.3  | 409 ±8     |
| C (HT. W/O CLOSURE):             | 9.60 ±.12   | 243.8 ±3.0 | 11.40 ±.14 | 289.6 ±3.6 | 12.90 ±.15 | 327.7 ±3.8 | 15.9 ±.2  | 404 ±5     |
| D (DIA).....                     | 4.75 ±.10   | 120.7 ±2.5 | 4.75 ±.10  | 120.7 ±2.5 | 6.00 ±.12  | 152.4 ±3.0 | 7.60 ±.13 | 193.0 ±3.3 |
| W.T.(NOM WALL THK):              | 0.07        | 1.8        | 0.08       | 2.0        | 0.10       | 2.5        | 0.12      | 3.0        |

**NOTES:**

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



**CONFIDENTIAL**

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher  
SCIENTIFIC**

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

|           |           |                    |           |
|-----------|-----------|--------------------|-----------|
| MADE MC   | APP'D PVC | PSS 00-220X-XXXX-1 | REV. 0403 |
| CHK'D PVC | 4/24/03   |                    |           |

DESCRIPTION: LARGE, NARROW-MOUTH ROUND BOTTLES