


CATALOG NO.	NOMINAL CAPACITY (GAL)	MATERIALS
002256	2 1/2	AMBER HDPE

SIZE CODE.....	-7020	
SIZE (NOM VOL L/GAL):	10/ 2.5	
NOM. BRIM VOL (L):	11.8	
NECK FINISH (mm):	83B	
THREAD (T.P.I.):	2	
NOM. WEIGHT W/ CLOSURE(gms):	1250	
	INCHES	MM
A (NECK I.D.):	2.52 ±.10	64.0 ±2.5
B (HT. W/ CLOSURE):	15.1 ±.20	383.5 ±5.0
C (HT. W/O CLOSURE):	14.7 ±.20	373.4 ±5.0
D (DIA):	9.85 ±.12	250.2 ±3.0
W.T.(NOM WALL THK):	0.12	3.0

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE MC ²	APP'D PVC	PSS 00-2256-7020	REV. 0105
CK'D PVC	1/5/05		

DESCRIPTION:
AMBER CARBOY WITH HANDLES