


CATALOG NO.	NOMINAL CAPACITY (L)	MATERIALS		
		CARBOY	CLOSURE	STRAP
002243-0013	13	HOPE	WHITE PPNAT	HOPE

SIZE (NOM VOL L).....	13	
NOM BRIM CAP (L).....	14.2	
NECK FINISH (mm).....	53	
THREAD (T.P.I.).....	4	
NOM. WEIGHT W/ CLOSURE(gms).....	1250	
	INCHES	MM
A (NECK I.D. AT TOP)..:	1.5	38
B (OVERALL HEIGHT)..:	14.9	378
C (WIDTH).....	11.4	290
D (DEPTH).....	7.45	189
W.T.(NOM WALL THK)..:	0.1	2.5



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE PB	APP'D PVC	PSS 00-2243-0013	REV. 1201
CK'D PVC	12/4/01		

DESCRIPTION:

13-LITER JERRICAN

2. This is not a production drawing, but is for informational purposes only, and is subject to change without notice.

1. All dimensions are nominal and in inches and millimeters unless otherwise noted.