


| CATALOG NO. | NOMINAL CAPACITY<br>(oz) |   |   |   |   | MATERIALS  |             |             |
|-------------|--------------------------|---|---|---|---|------------|-------------|-------------|
|             |                          |   |   |   |   | BOTTLE     | CLOSURE     | CAP         |
| 002411      | 1/2                      | 1 | 2 | 4 | 8 | LOPE       | PP          | PP          |
| 002414      | -                        | 1 | - | - | - | TEFLON FEP | TEFZEL ETFE | TEFZEL ETFE |

|                                  |           |           |           |           |           |
|----------------------------------|-----------|-----------|-----------|-----------|-----------|
| SIZE CODE.....:                  | -0015     | -0030     | -0060     | -0125     | -0250     |
| SIZE (NOM VOL mL/oz):            | 15/0.5    | 30/1      | 60/2      | 125/4     | 250/8     |
| NOM. BRIM CAP.(mL)::             | 18        | 35        | 64        | 140       | 285       |
| NECK FINISH (mm)::               | 20-415    | 20-415    | 20-415    | 24-415    | 24-415    |
| THREAD (T.P.I.)::                | 8         | 8         | 8         | 8         | 8         |
| NOM. WEIGHT<br>W/CLOSURE (gms):: | 8         | 10        | 12        | 17        | 28        |
| A (NECK I.D.)::                  | 0.54 ±.02 | 0.54 ±.02 | 0.54 ±.02 | 0.69 ±.02 | 0.69 ±.02 |
| C (HT. W/O CLOSURE)::            | 2.20 ±.06 | 2.30 ±.07 | 3.25 ±.08 | 3.88 ±.08 | 5.14 ±.09 |
| D (DIA)::                        | 0.98 ±.04 | 1.35 ±.05 | 1.52 ±.06 | 1.98 ±.06 | 2.40 ±.07 |
| H (NOM OVERALL HT.)::            | 3.8       | 3.9       | 3.9       | 4.5       | 4.8       |
| W.T.(NOM. WALL THK)::            | 0.040     | 0.040     | .035      | 0.040     | 0.050     |

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



NOT TO SCALE

SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Ponorona Creek Drive  
Rochester, NY 14625-2385, USA

|       |     |       |        |                    |      |      |
|-------|-----|-------|--------|--------------------|------|------|
| MADE  | MCC | APP'D | GRP    | PSS 00-241X-XXXX-1 | REV. | 0506 |
| CHK'D | DM  |       | 6/5/06 |                    |      |      |

DESCRIPTION:

DROP-DISPENSER BOTTLES