



CATALOG NO.	4110-0250	4110-0500		
NOM. BRIM CAP. (nL)	296	588		
MATERIAL	PC	PC		
NOM. WEIGHT	55	88		
	INCHES	MM	INCHES	MM
A (NECK I.D.)	1.06±.03	26.9±0.8	1.29±.02	32.8±0.5
B (NECK O.D.)	1.58±.03	40.1±0.8	1.49±.02	37.8±0.5
C (BASE DIA.)	3.25±.04	82.6±1.0	4.00±.04	101.6±1.0
H (HEIGHT)	5.26±.07	133.6±1.8	7.36±.07	186.9±1.8
W.T. (NOM. WALL THK)	0.05	1.3	0.05	1.3

CATALOG NO.	4110-1000	4110-2000		
NOM. BRIM CAP. (nL)	1120	2320		
MATERIAL	PC	PC		
NOM. WEIGHT	180	225		
	INCHES	MM	INCHES	MM
A (NECK I.D.)	1.20±.02	30.5±0.5	1.84±.03	46.7±0.8
B (NECK O.D.)	1.50±.02	38.1±0.5	2.27±.03	57.7±0.8
C (BASE DIA.)	5.11±.04	129.8±1.0	6.44±.05	163.6±1.3
H (HEIGHT)	9.64±.08	244.9±2.0	10.60±.08	269.2±2.0
W.T. (NOM. WALL THK)	0.06	1.5	0.06	1.5

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Pancromo Creek Drive
Rochester, NY 14625-2385, USA

NAME: MC
CR: B PVC

REV: D PVC 9/29/97

PSS 00-4110-XXXX

REV: 1200

DESCRIPTION:

BAFFLED FLASK