


CATALOG NO.	NOMINAL CAPACITY, (oz)						MATERIAL	
	1	2	4	8	16	16	32	BOTTLE
342024								PETG

SIZE CODE.....	-0030	-0060	-0125	-0250				
SIZE (NOM VOL nL/oz)	30/1	60/2	125/4	250/8				
NOM. BRIM CAP. (nL)	45	77	175	320				
NECK FINISH (mm)...	20-415	24-415	38-430	38-430				
THREAD (T.P.I.).....	8	8	6	6				
NOM. WEIGHT(gms)....	12.5	19.5	43.5	56				
	INCHES	MM	INCHES	MM	INCHES	MM	INCHES	MM
A (BTLE. NECK I.D.):	0.54±.02	13.7±.5	0.70±.02	17.8±.5	1.100±.025	27.94±.64	1.100±.025	27.94±.64
B (HT. BOTTLE):.....	2.39±.05	60.7±1.3	3.15±.06	80.0±1.5	4.08±.08	103.6±2.0	5.50±.08	139.7±2.0
C (SQUARE).....	1.50±.05	38.1±1.3	1.58±.05	40.1±1.3	2.06±.06	52.3±1.5	2.32±.06	58.9±1.5
MIN. WALL @ BTM COR.:	0.006	0.15	0.006	0.15	0.008	0.20	0.010	0.25
MIN. WALL ELSEWHERE.:	0.008	0.20	0.010	0.25	0.012	0.30	0.016	0.41

SIZE CODE.....	-0500	-9500	-1000			
SIZE (NOM VOL nL/oz)	500/16	500/16	1000/32			
NOM. BRIM CAP. (nL)	610	658	1225			
NECK FINISH (mm)...	38-430	38-430	38-430			
THREAD (T.P.I.).....	6	6	6			
NOM. WEIGHT(gms)....	107	107	178			
	INCHES	MM	INCHES	MM	INCHES	MM
A (BTLE. NECK I.D.):	1.100±.025	27.94±.64	1.100±.025	27.94±.64	1.100±.025	27.94±.64
B (HT. BOTTLE):.....	6.82±.09	173.2±2.3	6.82±.09	173.2±2.3	8.48±.10	215.4±2.5
C (SQUARE).....	2.88±.06	73.2±1.5	3.00±.06	76.2±1.5	3.63±.07	92.2±1.8
MIN. WALL @ BTM COR.:	0.014	0.36	0.014	0.36	0.014	0.36
MIN. WALL ELSEWHERE.:	0.020	0.51	0.017	0.43	0.020	0.51

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE TL	APP'D GRP	PSS 34-2024-XXXX	REV. 1009
CK'D DM 10/22/09	10/23/09		

PETG SQUARE MEDIA BOTTLES, STERILE