


CATALOG NO.	NOMINAL SIZE (mm)	MATERIAL CLOSURE	MATERIAL SEPTUM
342178	24	HDPE	SILICONE/PTFE

SIZE CODE	-0240	
THREAD FINISH (mm)	24-415	
THREAD (T.P.I.)	8	
NOM. WEIGHT (gms)	4	
	INCH	MM
D (DIA)	1.12 ±0.03	28.4 ±0.8
H (HEIGHT)	0.59 ±0.03	15.0 ±0.8
F (DIA)	0.47 ±0.03	11.9 ±0.8



SPECS ON BACK

NOT TO SCALE

**CONFIDENTIAL**  
 THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**  
 Thermo Fisher Scientific  
 75 Pennington Creek Drive  
 Rochester, NY 14625-2385, USA

MATERIAL	DATE	DESCRIPTION
MC <sup>2</sup> PVC	11/6/03	
APP'D PVC		

PSS 34-2178-0240

REV 1103

24-415 SEPTUM CLOSURE

**NOTES:**

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.