


CATALOG NO.....	3410-08	
SIZE NOM VOL. (L).....	10	
NOMINAL BRIM CAP. (L)....	13.6	
MATERIALS	BOTTLE.....	PC
	CLOSURE.....	WHITE PP
	CLOSURE LINER..	SILICONE RUBBER/PP
	HANDLE.....	LDPE
NECK FINISH.....	48mm DOUBLE START THR'D	
THREAD (T.P.I.).....	3	
NOM. WEIGHT W/ CLOSURE & HANDLE(gms):	770	
	INCH	CM
A (NECK I.D. MAX. TOP)...	1.49	3.8
A (NECK I.D. MIN. BTM)...	1.39	3.5
B (HT. W/O CLOSURE).....	13.17	33.5
C (HT. W/ CLOSURE).....	13.26	33.7
D (SQUARE).....	10.05	25.5
W.T. (NOM. WALL THK)....	0.08	0.2



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Ponorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MC <sup>2</sup>	APP'D	PVC	PSS 3410-08	REV.
CHK'D	PVC		1/9/06		0106

DESCRIPTION:  
BIOTAINER, CARBOY, PC, LAB, 10L WITH HANDLE