


CATALOG NO.	NOMINAL CAPACITY, (oz)						MATERIAL	
	1	2	4	8	16	16	32	BOTTLE
342024								PETG

SIZE CODE.....	-0030	-0060	-0125	-0250				
SIZE (NOM VOL nL/oz)	30/1	60/2	125/4	250/8				
NOM. BRIM CAP. (nL)	45	77	175	320				
NECK FINISH (mm)...	20-415	24-415	38-430	38-430				
THREAD (T.P.I.).....	8	8	6	6				
NOM. WEIGHT(gms)....	12.5	19.5	43.5	56				
	INCHES	MM	INCHES	MM	INCHES	MM	INCHES	MM
A (BTLE. NECK I.D.):	0.54±.02	13.7±.5	0.70±.02	17.8±.5	1.100±.025	27.94±.64	1.100±.025	27.94±.64
B (HT. BOTTLE):.....	2.39±.05	60.7±1.3	3.15±.06	80.0±1.5	4.08±.08	103.6±2.0	5.50±.08	139.7±2.0
C (SQUARE).....	1.50±.05	38.1±1.3	1.58±.05	40.1±1.3	2.06±.06	52.3±1.5	2.32±.06	58.9±1.5
MIN. WALL @ BTM COR.:	0.006	0.15	0.006	0.15	0.008	0.20	0.010	0.25
MIN. WALL ELSEWHERE.:	0.008	0.20	0.010	0.25	0.012	0.30	0.016	0.41

SIZE CODE.....	-0500	-9500	-1000			
SIZE (NOM VOL nL/oz)	500/16	500/16	1000/32			
NOM. BRIM CAP. (nL)	610	658	1225			
NECK FINISH (mm)...	38-430	38-430	38-430			
THREAD (T.P.I.).....	6	6	6			
NOM. WEIGHT(gms)....	107	107	178			
	INCHES	MM	INCHES	MM	INCHES	MM
A (BTLE. NECK I.D.):	1.100±.025	27.94±.64	1.100±.025	27.94±.64	1.100±.025	27.94±.64
B (HT. BOTTLE):.....	6.82±.09	173.2±2.3	6.82±.09	173.2±2.3	8.48±.10	215.4±2.5
C (SQUARE).....	2.88±.06	73.2±1.5	3.00±.06	76.2±1.5	3.63±.07	92.2±1.8
MIN. WALL @ BTM COR.:	0.014	0.36	0.014	0.36	0.014	0.36
MIN. WALL ELSEWHERE.:	0.020	0.51	0.017	0.43	0.020	0.51

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE TL	APP'D GRP	PSS 34-2024-XXXX	REV. 1009
CK'D DM 10/22/09	10/23/09		

PETG SQUARE MEDIA BOTTLES, STERILE