


PRODUCT DRAWING



- POLYESTER FILM WITH ACRYLIC ADHESIVE ON A CLEAR POLY-COATED LINER
- BACKSCORES THROUGH LINER ONLY
- PERFORATIONS THROUGH FILM ONLY. PERF = 0.2638 CUT/0.0182 TIE (6.7005-0.4623)

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. "( )" DIMENSIONS ARE IN MILLIMETERS. UNSPECIFIED TOLERANCES ARE ± 0.015.		DESC: ThermalSeal RT2RR™, TS-RT2RRNL-100			
DRAWN:	T.CARROLL	DATE:	11/14/13		
APPROVED:		DATE:	SCALE: 1:1	DWG #: MFG-DWG-410	REV:
			SHEET: 1 OF 1		