Edition Date: 2012-08-09 **Revision Date:** 2017-08-29

Revised by: Annette Roy, QA Coordinator

Product Name: EasyDip[™] Catalogue No.: M905-12DGY

This document replaces any previous version

1. Product Description:

➤ Slide Staining Rack: Non-sterile can hold up to 12 slides

> Compatible with:

• M900-12 EasyDip[™] Slide Staining System Series

• M906 Aluminum Holder

2. Packaging:

Case: 6 units per case (individually wrapped)

3. Product Specifications:

➤ Material:

Staining Rack: AcetalColor: Dark Gray

> Temperature range: 170°C to +80°C

4. Standards and Conformity:

> ISO 2859-1: Sampling and inspection procedures

CONEG / RoHS: Plastics and colorants are in conformity with

CONEG / RoHS standards for heavy metals

➤ **REACH (SVHC):** Plastic is in conformity to REACH standards

➤ LATEX: Material is Latex free

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension measurements

6. Traceability:

➤ Lot No. Composition: 8 or 9 digits

> The lot number can be found in one or all of these locations:

1. On exterior case label

2. On label inserted inside the master case

3. On the inner bag

Document Template: FRM-106 Rev01 Technical Data Sheet

Issue Date: 2017-06-15

Edition Date: 2012-08-09 **Revision Date:** 2017-08-29

Revised by: Annette Roy, QA Coordinator

7. **Storage Conditions:**

- > Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging

8. **Recommended Use:**

- > Avoid picric acid and phenol: prolonged use in these reagents will cause degradation
- > Follow manufacturer's instructions when using colorant or other chemical substance

OWNERSHIP OF MATERIALS:

Materials and information contained in this document are Simport Scientific's copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.

Page 2 / 2 Issue Date: 2017-06-15