


## Component Schematic

Krystal 24 Glass Bottom Microplate

P/N 324nnn



L	127.00 ± 0.25	X1	19.00 ± 0.10	C	6
W	85.65 ± 0.25	Xn	108.90 ± 0.10	R	4
H	20.00 ± 0.10	Y1	16.00 ± 0.10	P	18.00 ± 0.10
Hs	15.00 ± 0.10	Yn	70.00 ± 0.10	D	15.00 ± 0.10
Hw	18.45 ± 0.10	Hc	1.80 ± 0.10	after assembly	

Dimensions in mm

Well volume 3.1ml, after assembly

Well form Flat bottom

Material GPPS / Glass