

Safety at high speed

Improve the safety of your membrane filtration workflow

You don't have to be exposed to obvious risks every time you execute your microbiological analyses. Discover a safe and fast way of performing your membrane filtration with single-use EZ-Fit® Filtration Units.

Filtration set-up with reusable funnels


Flaming for disinfection can lead to burns or blistering.

Be careful when installing the funnel on the manifold, you might clip your fingers.




Flaming manifold head

5 sec



Flaming funnel base and inside

4 sec




Cooling time with sterile water

3 to 5 sec



Placing filter on manifold head

5 sec



Installing filter funnel on manifold

3 sec




Filtration set-up with sterile single-use EZ-Fit® Filtration Units



EZ-Fit® Filtration Unit set-up

1 sec



EZ-Fit® Filtration Units save you precious time while ensuring accurate results. The sterilized, disposable filtration devices for bioburden testing streamline your workflow.

Learn more:


[MerckMillipore.com/EZ-Fit](https://www.MerckMillipore.com/EZ-Fit)



For optimal microbial recovery, use EZ-Fit® Filtration Units with our Readyplate™ 55 plates. These 55 mm agar plates are small but mighty: the high filling volume ensures reproducible membrane transfer without bubbles and folds.

Learn more:

[MerckMillipore.com/ReadyPlate55](https://www.MerckMillipore.com/ReadyPlate55)



EZ-Fit® Filtration Unit, blue base with pad

Product Description	Pack Size	Cat. No.
EZ-Fit® Filtration Unit, white membrane, 0.45µm pore size, 100 mL funnel, bulk packaging	48	EFHAW100B
EZ-Fit® Filtration Unit, white membrane, 0.45µm pore size, 100 mL funnel, individual packaging	48	EFHAW100I
EZ-Fit® Filtration Unit, white membrane, 0.45µm pore size, 250 mL funnel, bulk packaging	48	EFHAW250B
EZ-Fit® Filtration Unit, white membrane, 0.45µm pore size, 250 mL funnel, individual packaging	48	EFHAW250I
EZ-Fit® Filtration Unit, black membrane, 0.45µm pore size, 100 mL funnel, bulk packaging	48	EFHAB100B
EZ-Fit® Filtration Unit, black membrane, 0.45µm pore size, 100 mL funnel, individual packaging	48	EFHAB100I
EZ-Fit® Filtration Unit, black membrane, 0.45µm pore size, 100mL funnel, bulk packaging	48	EFHAB250B
EZ-Fit® Filtration Unit, black membrane, 0.45µm pore size, 250 mL funnel, individual packaging	48	EFHAB250I



EZ-Fit® Filtration Unit, pink base

Product Description	Pack Size	Cat. No.
EZ-Fit® Filtration Unit, white membrane, 0.45µm pore size, 100 mL funnel, individual packaging, no pad	48	EFHVW10IS
EZ-Fit® Filtration Unit, white membrane, 0.45µm pore size, 100 mL funnel, multi packaging, no pad	48	EFHAW10MS
EZ-Fit® Filtration Unit, white membrane, 0.45µm pore size, 250 mL funnel, bulk packaging, no pad	48	EFHAW25BS
EZ-Fit® Filtration Unit, black membrane, 0.45µm pore size, 100 mL funnel, multi packaging, no pad	48	EFHAB10MS
EZ-Fit® Filtration Unit, black membrane, 0.45µm pore size, 250 mL funnel, bulk packaging, no pad	48	EFHAB25BS
EZ-Fit® Filtration Unit, white membrane, 0.22µm pore size, 100 mL funnel, multi packaging, no pad	48	EFGSW10MS
EZ-Fit® Filtration Unit, white membrane, 0.22µm pore size, 100 mL funnel, individual packaging, no pad	48	EFGSW10IS
EZ-Fit® Filtration Unit, white membrane, 0.8µm pore size, 100 mL funnel, bulk packaging, no pad	48	EFAAW10BS
EZ-Fit® Filtration Unit, white membrane, 0.8µm pore size, 250 mL funnel, bulk packaging, no pad	48	EFAAW25BS
EZ-Fit® Filtration Unit, black membrane, 0.8µm pore size, 100 mL funnel, bulk packaging, no pad	48	EFAAB10BS
EZ-Fit® Filtration Unit, black membrane, 0.8µm pore size, 250 mL funnel, bulk packaging, no pad	48	EFAAB25BS



EZ Family

Product Description	Pack Size	Cat. No.
EZ-Stream® Vacuum Pump for liquid transfer	1	EZSTREAM1
5 meters autoclavable silicone tubing with 9.5 mm internal diameter	1	STREAMTUB
EZ-Fit® Manifold, 1-place for EZ-Fit® Filtration Units	1	EZFITEFUN1
EZ-Fit® Manifold, 3-place for EZ-Fit® Filtration Units	1	EZFITEFUN3
EZ-Fit® Manifold, 6-place for EZ-Fit® Filtration Units	1	EZFITEFUN6
EZ-Fit® filtration head for EZ-Fit® Filtration Units	1	EZFITMVHE1
EZ-Fit® filtration head for EZ-Fit® Filtration Units	3	EZFITMVHE3
EZ-Fit® filtration head gasket	1	EZFITMVG3



ReadyPlate™ 55 agar plates

Product Description	Pack Size	Cat. No.
ReadyPlate™ 55 CCA ISO 9308	20x55mm	1467570020
ReadyPlate™ 55 CCA ISO 9308	20x55mm	1467570200
ReadyPlate™ 55 Lactose TTC Agar	20x55mm	1467590020
ReadyPlate™ 55 TGE Agar	20x55mm	1467610020
ReadyPlate™ 55 Plate Count Agar	20x55mm	1467630020
ReadyPlate™ 55 Slanetz and Bartley Agar ISO 7899	20x55mm	1467650020
ReadyPlate™ 55 CN Agar ISO 16266	20x55mm	1467670020
ReadyPlate™ 55 m-Green Agar	20x55mm	1467690020



To place an order or receive technical assistance in Europe, please call Customer Service:
 France: 0825 045 645 Spain: 901 516 645 Option 1
 Germany: 069 86798021 Switzerland: 0848 645 645
 Italy: 848 845 645 United Kingdom: 0870 900 4645

For other countries across Europe, please call: +44 (0) 115 943 0840
 Or visit: MerckMillipore.com/offices
 For Technical Service visit: MerckMillipore.com/techservice

© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, EZ-Fit, Millipore and the Vibrant M are trademarks of Merck KGaA, Darmstadt, Germany. All other trademarks are the property of their respective owners.

Lit. No. MK_FL2453EN 2018 - 13706 07/2018

